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ENV 50204: 1995, EN 61000-4-4: 1995, EN 61000-4-6: 1996)

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EN 61010

Thank you for purchasing a Heidolph Instruments product. This item has been designed, made and inspected in compliance with DIN EN ISO 61010 for long performance and continuous operation.

## SUMMARY

<b>SUMMARY .....</b>	<b>40</b>
<b>STANDARD ITEMS AND OPTIONS .....</b>	<b>42</b>
<b>GENERAL .....</b>	<b>43</b>
<b>SAFETY INFORMATION .....</b>	<b>43</b>
<b>SET-UP.....</b>	<b>45</b>
1. ROTAVAC valve control .....	45
1.1. Item set-up .....	45
1.2. Install cooler holder and condensate cooler .....	45
1.3. Vacuum supply .....	46
1.4. Electric branching .....	46
2. ROTAVAC vario control.....	46
2.1. Set-up .....	46
2.2. Install cooler holder and condensate cooler .....	46
2.3. Vacuum supply .....	47
2.4. Electric power supply .....	47
3. ROTAVAC vario tec and ROTAVAC valve tec .....	47
3.1. Set-up .....	47
3.2. Install cooler holder and condensate cooler .....	47
3.3. Vacuum supply .....	48
3.4. Electric power supply .....	48
4. VAC control automatic vacuum controller .....	49
4.1. Installation of VAC control automatic vacuum controller .....	49
4.2. Vacuum fitting .....	49
4.3. Power supply .....	49
5. Vacuum Switchbox .....	50
5.1. Set-up .....	50
5.2. Do connection work .....	50
<b>OPERATION AND CONTROLS.....</b>	<b>51</b>
6. ROTAVAC valve control, ROTAVAC valve tec, ROTAVAC vario control and ROTAVAC vario tec.....	51
6.1. Installation in rotary evaporator system: .....	51
6.2. Prior to first operation: .....	51
6.3. Operation: .....	51
6.4. Use gas ballast .....	52
6.5. Discontinue operation: .....	52
7. VAC control automatic vacuum controller .....	53
7.1. Configuration .....	53
7.2. Manual mode .....	54
7.3. Auto Mode .....	55
7.4. Temporary mode .....	56
7.5. Calibrate pressure sensor .....	57
7.5.1. Calibrate atmospheric pressure.....	57
7.5.2. Low value calibration .....	57

**INFORMATION ABOUT DISTILLATION SETTINGS ..... 58**

8. Rotary speed of evaporator flask ..... 58

9. Temperature difference between hot bath and vapor ..... 58

10. Temperature difference between vapor temperature and coolant ..... 58

11. Choose vacuum regulation mode and do respective settings ..... 58

11.1. "Manual" vacuum regulation mode ..... 58

11.1.1. Settings in "Manual" regulation mode ..... 58

11.2. "Auto Mode" vacuum regulation mode ..... 59

11.2.1. "Auto Mode" settings ..... 59

12. Solvent data ..... 59

12.1. at boil temperature of 40°C ..... 59

12.2. at boil temperatures deviating from 40°C ..... 59

12.3. solvents not depicted ..... 59

**CLEANING AND MAINTENANCE ..... 61**

13. ROTAVAC valve control ROTAVAC valve tec, ROTAVAC vario control and ROTAVAC vario tec ..... 61

13.1. Maintenance required ..... 61

13.2. Replace diaphragm and valve ..... 61

13.2.1. ROTAVAC valve control ..... 62

13.2.2. ROTAVAC vario control ..... 65

**DISASSEMBLE, FORWARD & STORE ..... 70**

**DISPOSAL ..... 70**

**TROUBLESHOOTING ..... 70**

**SPECIFICATIONS ..... 71**

**WARRANTY, LIABILITY & COPYRIGHT ..... 72**

**FAQ / REPAIR WORK ..... 73**

**CE-DECLARATION OF CONFORMITY ..... 73**

**Important information**

**Advice about power cord / mains supply**

**Caution: mandatory action**

**Caution: fire- and explosion hazard**

**Caution: hazard of burns, hot surface**

**Advice about maintenance / repair**

## STANDARD ITEMS AND OPTIONS

	item	qty.	P/N 230/240V 50/60Hz	P/N 115V 50/60Hz
	Rotavac valve control vacuum pump	1	591-00130-00	591-00130-01
or	Rotavac vario control vacuum pump	1	591-00140-00	591-00140-01
or	Rotavac valve tec vacuum pump		591-00160-00	591-00160-01
or	Rotavac vario tec vacuum pump		591-00170-00	591-00170-00
or	VAC control automatic vacuum controller	1	591-00340-00	591-00340-01
or	vacuum switchbox	1	591-00400-00	591-00400-01
<b>Standard items Rotavac valve control Rotavac valve tec</b>				
	AC power cord (valve tec)	1	14-007-002-78	14-007-002-79
	Owner's Manual	1	01-005-004-25	
	Warranty card	1	01-006-002-58	
<b>Standard items Rotavac vario control</b>				
	Rotavac vario control (power unit )	1	11-300-004-33	11-300-004-41
	Rotavac vario control (pump )	1	11-300-004-34	11-300-004-34
	Rotavac vario tec (pump )		11-300-004-73	11-300-004-73
	AC power cord	1	14-007-003-81	14-007-003-89
	Hose Connector Y-Type ID 8mm	1	23-30-01-03-52	
	Hose Nipple SO 20503-A8-8	2	23-09-01-01-79	
	Connection Cable 8 PIN	1	14-300-003-87	
	Owner's Manual	1	01-005-004-25	
	Warranty card	1	01-006-002-58	
<b>Standard items VAC control automatic</b>				
	Angle (Contr.autom.) lacquered	1	21-300-004-28	
	Nut M4	2	02-02-01-04-03	
	Washer A 4,3	6	02-03-06-04-04	
	Screw M4X12	4	02-01-04-02-04	
	Mains lead	1		14-300-003-86
	Hose Nipple PVDF	1	23-09-01-01-79	
	Hose Nipple M5	1	22-30-01-04-48-0	
		2	23-09-01-01-79	
	Owner's Manual	1	01-005-004-25	
	Warranty card	1	01-006-002-58	
<b>Standard items vaccum switchbox</b>				
	Non return valve	3	11-300-005-41	
	Hose connector Y-TYPE ID 8mmm	2	23-30-01-03-52	
	Connection cable (Vac senso switch)	3	14-300-003-67	
	AC power cord	1	14-007-003-81	14-007-003-89
	Owner's Manual	1	01-005-004-25	
	Warranty card	1	01-006-002-58	




### Options

	item	P/N
	Rotavac valve control condensate cooler ass'y.	591-00081-00
	Rotavac vario control condensate cooler ass'y.	591-00082-00
	Rotavac vario tec / valve tec condensate cooler ass'y.	591-00083-00

## GENERAL

-  **Unpack your item carefully.**  
Inspect for damage and report such damage or missing parts to your supplier immediately.
-  **Read your Instruction Manual carefully.** Take time to save time while working with your product. Make sure that every user has read and understood the Instruction Manual.
-  **Please store the Instruction Manual in a place easily accessible to every user.**

**IF ALL ELSE FAILS, READ THESE INSTRUCTIONS !**

-  A so-called EURO-plug (DIN 49441 CEE 7/VII 10/ 16 A 250 V) is standard on all of the products.  
For the Continental US they feature a US-standard plug (NEMA Pub.No.WDI.1961 ASA C 73.1 . 1961 page 8 15A 125V).
-  For using the item in a country with deviating outlet / plug systems, we recommend to use approved adapters or to have an electrician replace the standard plug with one mating your local system.
-  As shipped, the item features a protective ground wire. When replacing the original plug, make sure to reconnect this protective ground wire in the new plug !

## SAFETY INFORMATION

-  **Please comply with all safety and accident prevention regulations as in force for laboratory work !**
-  **Use extra care while working with vacuum pumps as used on rotary evaporators.**
-  **Use eye protection and adequate clothing.**
-  **Use extra care while working with flammable material. Refer to safety data sheets.**
-  **Use extra care while working close to flammable and explosive substances. Motors are of non-sparking type, the item itself is not explosion-protected.**

# E

Before connecting your item with the local power supply, make sure the item has been designed for your local voltage; refer to data plate on item.

Connect your item with a grounded outlet only.

Turn power switch to OFF when item is out of use, or before disconnecting from mains.

Refer to "hot surface" sticker on the pump.  
Do not touch hot surfaces.

Repair work is limited to skilled personnel so authorized by Heidolph.

Your item needs a proper stand !

Make sure no part of your body is exposed to vacuum pressure.

Avoid uncontrolled overpressure. Hazard of bursting ! Keep exhaust pipes clear of obstructions.

Due to non-compensated leak rate, a minor amount of gas exchange between ambient and vacuum system may occur. Preclude contamination of substances involved and environmental pollution.

At high intake pressures, due to high pump compression ratio, overpressure may occur at the gas ballast valve. Valve being open, gas or condensate may escape.

## SET-UP



## 1. ROTAVAC valve control

## 1.1. Item set-up

- Locate items on a plane, solid stand. Item must be horizontal.

## 1.2. Install cooler holder and condensate cooler

- Attach condensate cooler bracket (1) with 2 ea. screws from your hardware kit to the rear panel of the vacuum pump.
- Remove coupling nut (10) from condensate cooler (Attention, retainer underneath coupling nut!)
- Arrange condensate cooler (2) in holder.
- Secure condensate cooler (2) with rubber band (3)
- Install receiving flask (4) and secure with clamp (5).
- Remove exhaust stack (6) and replace with angular fitting (from hardware kit)
- Refer to drawing and install black coupling ring (8) on the left end of the PTFE hose (7); red coupling nut (10) and seal (9) go on the right-hand end of this hose. Refer to picture and slip this hose assembly over angular stud (6) and glass thread (18) on condensate cooler. Pull coupling rings (8 & 10) tight.
- Connect cooling water with fittings (11), exhaust air with fitting (14), as required.

## 1.3. Vacuum supply

- Connect vacuum with stud (12)

## 1.4. Electric branching

- Before connecting with your local power supply make sure that:

mains voltage is identical with data plate (refer to "safety Information"; refer to "General" if plug and outlet do not fit together.

- Connect power cord with outlet.

## 2. ROTAVAC vario control


**Caution:** ROTAVAC vario control may be used together with the LABOROTA 4002/4003 control vario evaporator only !

### 2.1. Set-up

- Locate items on solid and plain surface. Item must be horizontal.


### 2.2. Install cooler holder and condensate cooler

- If your ROTAVAC vario control is used together with the ROTACOOOL reflux cooler, you may directly place your ROTAVAC vario control on top of the ROTACOOOL.
- Using your ROTAVAC vario control without ROTACOOOL, an extra foot (15) has to be installed for the condensate cooler.
- Refer to picture and install extra foot (15) with 4 ea. screws (16), using hex pin from you hardware kit.
- Remove fairing (17) by unthreading two screws
- Install bracket (1) for condensate cooler (2) on the rear of the vacuum pump in a slot provided to this end. Introduce 2 ea. square nuts in slot.
- Re-install fairing (17).
- Remove coupling nut (10) from condensate cooler (Attention, retainer underneath coupling nut!)
- Locate condensate cooler (2) in holder
- Secure aftercooler with rubber band (3).
- Install receiving flask (4) and secure with clamp (5).
- Fit red coupling ring (10) and seal (9) on right-hand end of PTFE hose (7). Plug this line

assembly in glass thread (18) of aftercooler (2). Plug left-hand over outlet stud (6), pull coupling nut tight.

- Connect cooling water with fittings (11), exhaust air with fitting (14), as required.

### 2.3. Vacuum supply

- Vacuum supply is connected with inlet fitting (12)

### 2.4. Electric power supply

**Caution:** ROTAVAC vario control consists of a pump (24V=) and a power supply module (mains voltage)

- Before connecting power supply with mains network check whether:

actual voltage is identical with voltage depicted on data plate (refer to "Safety Information"); Refer to "General" for more information of plug connectors used.

- Connect power cord with wall outlet.
- Connect pump wire with control box (also refer to LABOROTA 4002/4003 control vario Owner's Manual)

## 3. ROTAVAC vario tec and ROTAVAC valve tec


**Caution:** ROTAVAC vario tec may be used together with the LABOROTA 4002/4003 evaporator only !

### 3.1. Set-up

- Locate items on solid and plain surface. Item must be horizontal.



### 3.2. Install cooler holder and condensate cooler

- If your ROTAVAC vario tec or ROTAVAC valve tes is used together with the ROTACOOOL reflux cooler, you may directly place your ROTAVAC vario tec or ROTAVAC valve tec on top of the ROTACOOOL.
- Remove fairing (17) by unthreading two screws
- Install bracket (1) for condensate cooler (2) on the rear of the vacuum pump in a slot provided to this end. Introduce 2 ea. square nuts in slot.
- Re-install fairing (17).

# E

- Remove coupling nut (10) from condensate cooler (Attention, retainer underneath coupling nut!)
- Locate condensate cooler (2) in holder
- Install receiving flask (4) and secure with clamp (5).
- Fit red coupling ring (10) and seal (9) on right-hand end of PTFE hose (7). Plug this line assembly in glass thread (18) of aftercooler (2). Plug left-hand over outlet stud (6), pull coupling nut tight.
- Connect cooling water with fittings (11), exhaust air with fitting (14), as required.

### 3.3. Vacuum supply

- Vacuum supply is connected with inlet fitting (12)

### 3.4. Electric power supply

- Before connecting power supply with mains network check whether:

actual voltage is identical with voltage depicted on data plate (refer to "Safety Information"); Refer to "General" for more information of plug connectors used.

**Caution:** ROTAVAC vario tec consists of a pump (24V=) and a power supply module (mains voltage)

- Connect power cord with wall outlet.
- Connect pump wire with control box (also refer to LABOROTA 4002/4003 Owner's Manual)


## 4. VAC control automatic vacuum controller

### 4.1. Installation of VAC control automatic vacuum controller

The VAC control automatic vacuum controller is installed on top of the lift column of the rotary evaporator. Item requires little space, control panel is accessible at ease.

- refer to picture and install metal bracket (20) on vacuum controller 22; use two ea. nuts (21).
- install assembly on the drive head's (24) rear panel with 4 ea. screws (23) and hex. pin (25).


### 4.2. Vacuum fitting

refer to picture

### 4.3. Power supply

refer to picture

- A power switch
- B vacuum valve
- C mains
- D vacuum
- E vacuum switchbox
- F aeration


## 5. Vacuum Switchbox

### 5.1. Set-up

Locate vacuum switchbox on a solid and plain surface.

### 5.2. Do connection work

This vacuum switchbox allows to run up to 3 ea. rotary evaporators with one only vacuum pump. For branching and vacuum supply refer to picture. Since these 3 ea. vacuum evaporators may be run at different pressures, installation of 3 ea. non-return valves is mandatory. Lines and non-return valves come with the item. Connect vacuum pump with plug connector of the vacuum switchbox.

## OPERATION AND CONTROLS

Before connecting the item with your local power supply make sure that:

- the item has been designed for your local power supply network (refer to data plate on the item)

## 6. ROTAVAC valve control, ROTAVAC valve tec, ROTAVAC vario control and ROTAVAC vario tec

### 6.1. Installation in rotary evaporator system:

- Avoid obstruction losses by using short vacuum lines of high cross section.
- Avoid transfer of mechanical loads by stiff lines; use elastic hoses or shock attenuators instead. Caution: elastic elements may contract during evacuation.
- Avoid condensate reflux from the pump's exhaust end; exhaust pipes need to be routed with adequate slope.

### 6.2. Prior to first operation:

- Max. ambient temperature: 40 °C.
- Installing in whatever sort of housing or while working at high ambient temperatures make sure to have adequate ventilation. Maintain a safety offset of 20 cm between fans and adjacent parts.
- Gas ballast valve of the pump being open, electric power failure may cause aeration. If such situation might be hazardous, take adequate action.
- Avoid whatever sort of leakage in your installation; check system for leaks after installation.

### 6.3. Operation:

The pump may be started at 1 bar max. differential pressure between in- and outlet port to avoid damage to the motor.

- Refer to max. admissible pressure at outlet port and max. admissible pressure difference between in- and outlet port of the pump (refer to "Specifications" section).

High compression ratio of the pumps may cause high pressure at output end; never close the output end of the pump !

The pump is going to reach specified suction and output pressure values, as well as vapor compatibility at normal operating temperature (after abt. 15 minutes).

- Avoid condensation in the pump, hard fluid pulses and dust; continuous feeding of fluids or dust may damage diaphragm and valves.
- Run pump with gas ballast to minimize condensation of pumped substances in the pump (steam, solvents, .....).

A thermal winding protector will disconnect motor power supply in case of overtemperature.


- Caution: accepts manual reset only. Turn pump OFF and disconnect power plug. Allow about 5 minutes for cooling before restart.

## 6.4. Use gas ballast

- Make sure that air / gas entry through the ballast valve does not produce reactive, explosive or otherwise hazardous mixes. In case of doubt use inert gasses. If not sure, get in touch with Heidolph.
- Using air instead of inert gasses may involve damage to your hardware and/or the environment, serious injuries, danger to life due to production of hazardous and / or explosive mixes, if air and substance pumped react in the pump or at the pump outlet port.

Condensable vapors (steam, solvents, ...):

- Open gas ballast valve (refer to picture). Rotate through 180° to close valve.
- The gas ballast valve being open, final vacuum specified will not be reached.
- Never load pump with steam before having reached operating temperature.

Condensate produced:

- Use Woulff flask (optional) and condensate cooler (optional) (no condensate reflux, controlled condensate collecting).

Not using a condensate cooler at the output end, hazardous substances may escape; refer to your local regulations !

Check fluid level in Woulff flask and receiving flask of the condensate cooler; empty on condition:

- open clamp, remove condensate cooler, drain condensate.
- aerate system via vacuum fitting, unthread all hoses, and remove Woulff flask, drain condensate.

Caution: dispose of chemical substance in accordance with local regulations.

## 6.5. Discontinue operation:

### Short-term:

Is there any accumulation of condensate in the pump ?

Allow the pump to run a few extra minutes with open intake stud.

Is there any accumulation of substances in the pump that might degrade pump material ?

- Clean and check pump heads on condition.


### Long term:

- Do actions as described for short-term putting out of operation.
- Disconnect pump from system.
- Close manual gas ballast valve.
- Close inlet and outlet ports (e.g. with caps).
- Store pump in a dry place.

## 7. VAC control automatic vacuum controller

Power switch is located on the right.

When turning **VAC control automatic** ON, the right-hand display (Act Vac) shows the unit of

pressure (hPa=mbar) for a few seconds. Thereafter, modes, nominal value and actual value are displayed (both in mbar).

Start control action with the "Auto" or "Man" keys.

### 7.1. Configuration

The VAC control automatic comes with the following configuration:

pressure values = mbar (hPa)

manual mode = H1

automatic-parameters = "35"

Edit this basic configuration:

turn item OFF; depress "Auto" key and hold depressed while turning the item OB again.

the "Act Vac" display first shows the unit chosen. Using the arrow keys, you may either select "hPa" (mbar) or "Torr". We recommend to use the "hPa" (mbar) mode.

depress "Set" key, "H1" is displayed. You may change to "H2" with the arrow keys. "H1" is the standard regulating mode. vacuum valve is opened and closed normally. "H2" is going to open the vacuum valve for a second only; this way it is easier to maintain nominal values.

Using a switchbox, "H1" is the only operational mode available.

# E

depress "Set" key, "A1" is displayed; at the right of "A1" a numeric value for the automatic parameter appears, e.g. 35 (refer to Chapter 6.3). This value may be changed with the arrow keys between 20 (sensitive) and 50 (inert).


depress "Set" key, configuration is complete.

you may quit the configuration mode at whatever point with the "STOP" key.

## 7.2. Manual mode

In manual mode, a manually defined nominal pressure is used. The following nominal profile may be selected:

- 1: rapid decrease to nominal
- 2: ramp-shaped pressure decrease (refer to differential pressure and decrease interval)
- 3: end of cycle, vacuum valve closed, aeration valve open.

The following parameters may be set:

- nominal pressure, 1-999 mbar
- hysteresis (h), 1-100 mbar
- differential pressure, decrease:  $\Delta p$  (d), 0-100 mbar
- decrease time:  $\Delta t$  (t), toff (without) – 99 min

The "Set Vac" display is used for dialog display. All values may be edited with the arrow keys. After not having done any editing for ten seconds, the item returns to "Actual".

Example:

The "Act Vac" display shows the nominal value to be set (display flashing); for setting, use arrow keys.

The "Act Vac" display shows the hysteresis to be set, e.g. h\_\_5; this value may be adjusted with the arrow keys.

The "Act Vac" display shows differential pressure to be set, e.g. d\_45 which in turn may be edited with the arrow keys; no differential pressure wanted, enter "0".

The "Act Vac" display shows the time decrease to be set, e.g. t\_15, which in turn may be edited with the arrow keys. Not needing the "decrease pressure" mode, select „t off“ with the arrow keys.

Return to actual value display.

Start and discontinue cycle:

start regulation, mode display shows "H1" (LED "MAN" lighted)

Running a decrease ramp, time remaining will alternate with the "H1" mode display. At cycle end, "END" will be displayed, valves are closed, and a buzzer sounds.

Depressing the "Man" key a second time during cycle, the item will run against final vacuum. "PU" is displayed.

Discontinue regulation, aerate rotary evaporator.

### 7.3. Auto Mode

In the automatic mode, first boil point ("bP": boiling point) is found automatically while decreasing pressure.

For optimization, an automatic parameter is available (refer to Chapter 7.1). A low value selected will go in hand with a high-sensitive automatic reaction. A value being too low, an early boil point reaction may be the result. A value being too high may cause "overriding" boil point.

As in manual mode, a pressure decrease ramp may be entered (refer to Chapter 6.1), which is used automatically after having detected the first boil point.

The following parameters may be selected:

- nominal (unimportant)
- hysteresis (h), 1-100 mbar
- differential pressure, decrease (d), 0-100 mbar
- decrease time (t), toff (without) – 99 min

Enter parameters as described in Chapter 7.1.

Start and discontinue cycle:

start regulation, "A1" appearing in mode display, AUTO LED lighted.

Arrow and SET keys are locked; after having found the first boil point, "bP" is displayed.

Having programmed a decrease ramp, time remaining will alternate with "bP" display. At cycle end, "End" is displayed, and all valves are closed.

No decrease ramp having been programmed, boil pressure once determined will be controlled until the cycle is discontinued with "STOP".

Discontinue regulation, aerate rotary evaporator.

After having found a first boil point, you may change to manual operation whenever you want to do so. To this end:

depress "Man" key. LED beside "Man" key is flashing.

depressing the "SET" key, "bP" boil point will be saved as new nominal value.

## 7.4. Temporary mode

Independent of program cycle, the user may step into such cycle whenever he wants ("Manual" mode), and set a temporary nominal value, which will not be saved as operational parameter.:

by briefly hitting an arrow key, the present actual value will be set as temporary actual value ("MAN" LED flashing). Such temporary nominal value may be edited with the arrow keys.

Holding this key depressed, the system will be evacuated at full pumping power. Nominal value continuously is governed by present actual pressure.

Holding this key depressed, the aerating valve is opened 3 times, followed by continuous opening. Nominal value continuously is governed by present actual pressure.

Exit temporary mode with "Man" key; original nominal value is recalled.

Discontinue regulation, aerate rotary evaporator.

## 7.5. Calibrate pressure sensor

The VAC control automatic pressure sensor may be calibrated.

To this end you need an accurate vacuum gauge which is installed in the vacuum line between vacuum pump and VAC control automatic vacuum fitting.

### 7.5.1. Calibrate atmospheric pressure

Remove vacuum hose from VAC control automatic.

Turn VAC control automatic power switch OFF. Depress "Set" key, hold depressed while turning power switch ON. The "Act Vac" display shows your firmware's version number.

Depress "Set" key, "CAL" is displayed.

Depress arrow key "Up", present atmospheric pressure appears in the "Act Avp" display; "Set Vac" display shows "HI"

Copy pressure gauge reading with arrow keys

Complete calibration by hitting the "Set" key.

### 7.5.2. Low value calibration

Connect vacuum hose with VAC control automatic and start vacuum pump. Pressure decreases to lowest possible level as achieved by the vacuum pump.

If vacuum pump does not reach a pressure of less than 100 mbar, calibration is discontinued; an error report appears in the Act Vac display.

Turn power switch OFF. Depress "Set" key, hold and turn power switch ON. The "Act Vac" display shows your firmware's version number.

Depress "Set" key, "CAL" is displayed.

Depress arrow key "down", present pressure appears in the "Act Vac" display. "Set Vac" display shows "LO".

# E

Copy vacuum gauge reading with arrow keys.

Complete calibration by hitting the "Set" key

## INFORMATION ABOUT DISTILLATION SETTINGS

Using a rotary evaporator for thermal cracking of substances, some basic considerations should be made to achieve the best results possible. Optimize settings as described below:

### 8. Rotary speed of evaporator flask

Raise distillation speed by raising rotational speed. This method cuts distillation time and protects the substances handled.

### 9. Temperature difference between hot bath and vapor

This temperature difference should at least be 20 K to achieve adequate distillation speed. If acceptable by thermal stability of the substances handled, you even may raise temperature difference.

Rule of thumb: double temperature to double distillation speed.

### 10. Temperature difference between vapor temperature and coolant

This temperature difference should be no less than 20 K to guarantee adequate condensation. With low temperature difference, effectiveness of solvent recovery drops. Choose an as high as possible temperature difference in particular for high-boiling substances. .

### 11. Choose vacuum regulation mode and do respective settings

#### 11.1. "Manual" vacuum regulation mode

The "Manual" mode features the following parameters:

- separation of substances in a mix within physical limits
- no optional equipment needed
- manually select pressure when distilling mixes of substances
- high repetition accuracy
- requires detailed knowledge about thermal dynamic substance parameters

particularly for all standard processes

#### 11.1.1. Settings in "Manual" regulation mode

- pressure "p" (Set Vac) and the hysteresis  $\Delta p$  (Set  $\Delta p$ ) are set.
- choose pressure "p" as to locate boiling point of the substance to be distilled in the limits defined by Chapters 9 and 10 (bath temperature and coolant temperature). To this end make

reference to reference chart (page 60) and nomogram (page 60).

- hysteresis  $\Delta p$  defines temperature rise in mbar, required to open the vacuum valve.

Extended valve opening times equals high losses due to repeated gas transport.

- Low-boiling substances prefer a higher  $\Delta p$  (5 - 10 mbar), whereas we recommend a lower value (1 - 5 mbar) for high-boiling substances.

## 11.2. "Auto Mode" vacuum regulation mode

The "Auto Mode" is distinguished by the following parameters:

- Find distillation point automatically
- Distillation pressure once found will be maintained permanently, integrating hysteresis  $\Delta p$  (h).
- Separate mixes within physical limits
- No need to know about thermal dynamic substance data.

### 11.2.1. "Auto Mode" settings

- The only item to be set is h hysteresis  $\Delta p$  (h).

## 12. Solvent data

Examples for using reference table and nomogram:

This reference table details substance data needed for distillation, whereas the nomogram reflects relationship between pressure and boil temperature of some solvents. Go by the Clausius-Clapeyron equation and plot  $1/T$  against  $\log p$ .

### 12.1. at boil temperature of 40°C

LABOROTA 4002/4003 control: in column "vacuum for boil temperature" use the "Set Vac" nominal value as depicted for 40°C.

### 12.2. at boil temperatures deviating from 40°C


4. On temperature axis of the nomogram mark boil point wanted.
5. Draw a line to the right and determine intersection with the solvent line.
6. From this intersection draw a straight line down to determine the vacuum pressure needed.

### 12.3. solvents not depicted

- The following items might help in determining the right vacuum pressure:
  1. Rise of the straight line is a function of evaporation enthalpy. Hence, straight lines depicted may be used for rough reference to a lightly deviating boil point.
  2. A water jet pump or diaphragm-type pump may lower boil temperature by about 100°C.
  3. Rule of thumb: decrease pressure to  $\frac{1}{2}$  will decrease boil point by about 15 %.

Solvent	Total formula	MW [g/mol]	MW [°C]	$\Delta H_{vap}$ [J/g]	Vacuum for Boiling point at 40°C	
					[mbar]	[mm(Hg)]
Aceton	C <sub>3</sub> H <sub>6</sub> O	58,08	56,5	550	556	387
Acetonitril	C <sub>2</sub> H <sub>3</sub> N	41,05	81,8	833	230	173
Benzol	C <sub>6</sub> H <sub>6</sub>	78,11	80,1	549	236	177
n-Butanol (Butylalkohol)	C <sub>4</sub> H <sub>10</sub> O	74,12	117,5	619	25	19
tert.-Butanol (tert.-Butylalkohol)	C <sub>4</sub> H <sub>10</sub> O	74,12	82,9	588	130	98
2-Butanon (Methylethylketon)	C <sub>4</sub> H <sub>8</sub> O	72,11	79,6	473	243	182
tert.-Butylmethylether	C <sub>5</sub> H <sub>12</sub> O	88,15	55,0			
Chlorbenzol	C <sub>6</sub> H <sub>5</sub> Cl	112,60	132,2	375	36	27
Cyclohexan	C <sub>6</sub> H <sub>12</sub>	84,16	80,7	389	235	176
1,2-Dichlorethan	C <sub>2</sub> H <sub>4</sub> Cl <sub>2</sub>	98,96	82,4	336	210	158
1,2-Dichlorethylen (cis)	C <sub>2</sub> H <sub>2</sub> Cl <sub>2</sub>	96,94	59,0	320	479	134
1,2-Dichlorethylen (trans)	C <sub>2</sub> H <sub>2</sub> Cl <sub>2</sub>	96,94	47,8	313	751	563
Dichlormethan (Methylenchlorid)	CH <sub>2</sub> Cl <sub>2</sub>	84,93	40,7	373	atm.	atm.
Diethylether	C <sub>4</sub> H <sub>10</sub> O	74,12	34,6	392	atm.	atm.
Diisopropylether	C <sub>6</sub> H <sub>14</sub> O	102,20	67,5	318	375	281
Dimethylformamid	C <sub>3</sub> H <sub>7</sub> NO	73,09	153,0		11	8
1,4-Dioxan	C <sub>4</sub> H <sub>8</sub> O <sub>2</sub>	88,11	101,1	406	107	80
Ethanol	C <sub>2</sub> H <sub>6</sub> O	46,07	78,4	879	175	131
Ethylacetat	C <sub>4</sub> H <sub>8</sub> O <sub>2</sub>	88,11	77,1	394	240	180
Heptan	C <sub>7</sub> H <sub>16</sub>	85,09	98,4	439	120	90
Hexan	C <sub>6</sub> H <sub>14</sub>	86,18	68,7	370	335	251
Methanol	CH <sub>4</sub> O	32,04	64,7	1225	337	253
3-Methyl-1-Butanol (Isoamylalkohol)	C <sub>5</sub> H <sub>12</sub> O	88,15	130,6	593	14	11
Pentachlorethan	C <sub>2</sub> HCl <sub>5</sub>	202,30	160,5	203	13	10
Pentan	C <sub>5</sub> H <sub>12</sub>	72,15	36,1	382	atm.	atm.
n-Pentanol (Amylalkohol)	C <sub>5</sub> H <sub>12</sub> O	88,15	137,8	593	11	8
1-Propanol (n-Propylalkohol)	C <sub>3</sub> H <sub>8</sub> O	60,10	97,8	787	67	50
2-Propanol (Isopropylalkohol)	C <sub>3</sub> H <sub>8</sub> O	60,10	82,5	701	137	103
1,1,2,2-Tetrachlorethan	C <sub>2</sub> H <sub>2</sub> Cl <sub>4</sub>	167,90	145,9	247	35	26
Tetrachlorethylen	C <sub>2</sub> Cl <sub>4</sub>	165,80	120,8	233	53	40
Tetrachlormethan (Carbontetrachlorid)	CCl <sub>4</sub>	153,80	76,7	225	271	203
Tetrahydrofuran	C <sub>4</sub> H <sub>8</sub> O	72,11	66,0		357	268
Toluol	C <sub>7</sub> H <sub>8</sub>	92,14	110,6	425	77	58
1,1,1-Trichlorethan	C <sub>2</sub> H <sub>3</sub> Cl <sub>3</sub>	133,40	74,1	251	300	225
Trichlorethylen	C <sub>2</sub> HCl <sub>3</sub>	131,40	86,7	265	183	137
Trichlormethan (Chloroform)	CHCl <sub>3</sub>	119,40	61,3	263	474	356
Wasser	H <sub>2</sub> O	18,02	100,0	2259	72	54
Xylol (Isomeren-Gemisch)	C <sub>8</sub> H <sub>10</sub>	106,20	137-143	390	25	19

Conversion factor [mbar] to [mm(Hg)]:  $[mm(Hg)] \approx \frac{3}{4} [mbar]$

## CLEANING AND MAINTENANCE

For cleaning, wipe housing with a damp cloth (add some sort of mild liquid soap)..

**NOTE:** To avoid damage to surface finish, avoid using chlorine bleach, chlorine-based cleaners, abrasive substances, ammonia, rags or cleaning pads containing metallic particles.

VAC control automatic and Switchbox run maintenance free. Repair work is limited to engineers authorized to do so by HEIDOLPH. To this end get in touch with your local HEIDOLPH dealer or field representative.

### 13. ROTAVAC valve control ROTAVAC valve tec, ROTAVAC vario control and ROTAVAC vario tec

#### 13.1. Maintenance required

All Heidolph Vacuum Pumps require cleaning and maintenance to clean out chemical build up from solvents not distilled. We recommend a weekly rinse, or more depending on use, with Acetone or equivalent to break up and clean out solvent.

Disconnect all tubes and switch ON the pump. Slowly drip Acetone into the INLET nozzle of the operating pump. The cleaning solvent will be distributed through the pump valves and along the diaphragms, removing remains. Please ensure that your cleaning solvent exiting the pump is collected in a vessel at the OUTLET nozzle.

**ATTENTION:**

Be careful not to use too large of a quantity of cleaning solvent in the pump, which will cause damage to the pump diaphragms.

#### 13.2. Replace diaphragm and valve

All bearings are sealed hermetically and do not need any more lubrication. The pump itself is maintenance-free. Valves and diaphragms are wear parts. As soon as pressure drops you should clean and check housing, diaphragms and valves for cracks.

From case to case it make sense to clean and inspect pump heads in regular intervals. Normal service life of diaphragms and valves is > 8000 hours.

- Continuous feeding of fluids and dust will harm diaphragm and valves. Avoid condensation in pump, as well as accumulation of dust.

Pumping corrosive gasses and vapors or in case of deposits in the pump, such maintenance work should be done more frequently (by user's experience).

- Regular maintenance will increase pump's service life and protect operator and environment.

Before connecting maintenance work, disconnect pump from system, and pull mains plug. Allow pump to chill. Take adequate action to avoid pump running in disassembled condition.

Did you pump hazardous or corrosive gasses ?

- Take precautions (protective garments and goggles), to avoid breathing and skin contact.

### 13.2.1 ROTAVAC valve control

Seal kit (diaphragms and valves) for ROTAVAC valve control valve  
diaphragm  
special spanner, size 40/4

11-300-005-39  
23-30-01-01-89  
23-30-01-01-90  
02-07-02-01-12


#### Tooling:

- Phillips screwdriver # 2
- open end spanner, SW 15/17
- hex pin, SW 5
- special spanner, size 40/4


#### Clean and inspect pump heads:

Loosen coupling rings with open end spanner (SW 17).


Remove hose fitting from hose with box-end spanner (SW 15) by  $\frac{1}{4}$  revolution. Do not unthread fittings from the pump heads.

- Check for leaks after re-installation.


**Explosion drawing: pump head**


Loosen 4 ea. machine hex-recessed bolts on pump heads and lift off upper half (lid with internal part)

- Never attempt to loosen parts with sharp tools (screwdriver), use rubber hammer or compressed air instead.


To inspect valves carefully detach lid from internal part; denote location of valves and remove valves.

- Replace damaged valves.
- Remove contamination with cleaning agent.  
 Inspect diaphragm for damage, replace on condition. To this end remove data plate with 4 ea. Phillips screws.

## Replace diaphragm:

Loosen diaphragm back-up disc with special tool.

Check for shims; arrange shims by cylinder and re-install same number of shims.

- Insufficient number of shims: pump won't be able to reach final pressure; too many shims: pump running noisy.

Install new diaphragm between back-up disc with square-head bolt and back-up ring.

- Caution: light-colored side of diaphragm points towards cocking disc.


Install cocking disc, diaphragm, and back-up ring with piston rod; use special tool.

- Make sure square-head bolt is properly seated in the cut-out of back-up ring.


## Install pump heads:

Turn eccentric until piston rod reaches a position that allows diaphragm to be located in the center of the diaphragm and smoothly on mating surface of the housing's cut-out.

Reverse above sequence of operations for assembly. Install cover, valves and inside part of lid and cover.

- Make sure to install valves the right way.
- Turn eccentric until piston rod reaches upper dead center.
- Max. piston rod stroke.
  - Pull 4 ea. hex-recessed straight screws diagonally (hand-tight) (e.g. in the following sequence ←, ↑, →, ↓), thereafter torque to 10 Nm.


## Check serviceability of pump heads

Gauge pressure at intake port.

With an adequate vacuum gauge reading at the intake end should be less than 90 mbar.


- If not, check housing once more. Check correct seating of valves and diaphragms (in-line with hole in housing).


**Install flexible connection:**

turn fitting into hose with box-end spanner (SW 15).

Pull tight coupling ring by hand, thereafter pull tight with one revolution (use special tool).

**Pump does not reach final pressure specified:**

- After having replaced diaphragm or valves, pump will need several hours until reaching the pressure specified.

If after having replaced diaphragm or valves, actual values are off nominal, nothing will change by bleed-in.

First inspect screwed fitting of hoses on pump heads. As required, loosen screwed fitting and do one more revolution on elbow.

**13.2.2. ROTAVAC vario control**

Seal kit ROTAVAC vario control (4 ea. diaphragms, 8 ea. valves, special spanner) 11-300-005-40

Seal kit ROTAVAC vario tec and ROTAVAC valve tec (2 ea. diaphragms, 4 ea. valves, special spanner) 11-300-009-16

valve 23-30-01-01-91

diaphragm 23-30-01-01-92


special spanner SW 46 02-07-02-01-13

**tooling:**

- special spanner SW 46 (comes with seal kit)
- box end spanner SW 15/17
- hex pin SW 4/5
- Phillips screwdriver size #2

**Explosion view: pump head**

**Clean and inspect pump heads**

**Only ROTAVAC vario control:**

Use box-end spanner (SW 17) and loosen coupling ring of hose beside gas-ballast valve. Remove hose fitting from hose with ¼ revolution with box-end spanner (SW 15). Do not remove fittings from pump head.

- Leaks may be found on re-assembly.

**ROTAVAC vario control, ROTAVAC vario tec und ROTAVAC valve tec:**


Loosen 6 ea. straight hex-recessed screws on pump head and lift off cover ass'y. (cover, internal part, and lid):

- Never attempt to loosen parts with sharp tools (screwdriver), use rubber hammer or compressed air instead.

To inspect valves, carefully remove cover from lid. Denote location of valves and remove.

- Replace damaged valves.
- Remove contamination with cleaning agent.  
Inspect diaphragm for damage and replace on condition.

## Replace diaphragm

Lift diaphragm carefully.


- As required, move piston rod to upper dead end by depressing the cocking disc.
- Do not use sharp tools to lift diaphragm.

Introduce special tool underneath diaphragm.

Depress cocking disc until reaching lower dead center.


Force special tool against cocking disc and unscrew cocking disc with diaphragm.

- If old diaphragm wont detach from back-up disc, try to us gasoline or petroleum.
- Check for shims originally installed, use the same number for reassembly.

Introduce new diaphragm between cocking disc and back-up disc.

- Caution: Light-colored side of diaphragm points towards cocking disc.


Lift edge of diaphragm and carefully introduce cocking disc and back-up ring in special tool.

- Avoid damage of diagram: do not kink diaphragm until showing light-colored lines on the surface.

Check for shims originally installed, do not remove, re-install same number.

- Insufficient number of shims: pump won't reach final pressure; too many shims: pump running noisy.


# E

## Install pump heads

- Correct location of the square-head screw of the cocking disc in the slot of the back-up disc is of importance.


Screw diaphragm, cocking disc and back-up ring with piston rod.

- Insert shims between back-up ring and piston rod as required.


Locate diaphragms in a central and plain position on the housings cut-out.

Lay down pump and support as required.


Install cover and valves.

- Check for proper location (refer to fig.)

Proper location and orientation of cover and valves are mandatory.


Install lid.

- Make sure lid is seated properly.
- Pull 6 ea. hex. recessed bolts hand-tight diagonally,
- thereafter torque to 6 Nm.

**Install flexible connection (only ROTAVAC vario control):**

Install fitting in hose with box-end spanner (SW 15).  
Pull coupling ring right by hand; thereafter make one more turn with box-end spanner.

**Pump does not reach final pressure specified:**

- After having replaced diaphragm or valves, pump will need several hours until reaching the pressure specified.
- If after having replaced diaphragm or valves, actual values are off nominal, nothing will change by bleed-in.

# E

Inspect screwed fitting of hoses on pump heads.

## DISASSEMBLE, FORWARD & STORE

### Disassemble

Turn item OFF and disconnect mains plug.

### Transportation and store

1. Store item and its components in the original box to avoid damage in transit. Close box with adhesive tape.
2. Store item in a dry place.

**Caution: avoid shocks and vibration in transit.**

## DISPOSAL

For disposal, please comply with your local or national regulations.  
 Split by metal, plastic, etc.  
 Packing material to be treated as described above (material split).

## TROUBLESHOOTING

### Vacuum pumps

defect found	reason	correction
pump won't run or starts running and stops right away	<ul style="list-style-type: none"> <li>➤ electricity missing, voltage too low ?</li> <li>➤ <b>ROTAVAC vario control and ROTAVAC vario tec:</b> triggering signal for motorspeed missing?</li> <li>➤ overpressure in exhaust line?</li> <li>➤ motor overload?</li> </ul>	<ul style="list-style-type: none"> <li>➤ check voltage, connect electricity</li> <li>➤ check triggering signal</li> <li>➤ open exhaust line.</li> <li>➤ allow motor to chill, check for reason, turn ON again</li> </ul>
low suction power	<ul style="list-style-type: none"> <li>➤ centering ring not seated correctly</li> <li>➤ long, thin line?</li> <li>➤ condensate in pump?</li> <li>➤</li> <li>➤ deposits in pump ?</li> <li>➤ diaphragm and valves broken?</li> <li>➤ emission of gas, vapor in process ?</li> <li>➤ <b>ROTAVAC vario control and Rotavac vario tec:</b></li> <li>➤ pump runs hot (decrease speed)?</li> </ul>	<ul style="list-style-type: none"> <li>➤ check pump check fittings. line and recipient</li> <li>➤ increase cross section.</li> <li>➤ run pump with intake stud open for a minute</li> <li>➤ clean and check pump heads</li> <li>➤ replace diaphragm and/or valves</li> <li>➤ check process parameters</li> <li>➤ improve cooling, or reduce suction pressure</li> </ul>

pump noisy	<ul style="list-style-type: none"> <li>➤ sucks at atmospheric pressure ?</li> <li>➤ cocking disc loose?</li> <li>➤ none of the above reasons found ?</li> </ul>	<ul style="list-style-type: none"> <li>➤ install hose on exhaust end</li> <li>➤ maintain diaphragm pump</li> <li>➤ ship pump for repair.</li> </ul>
pump stuck or piston rod jammed		- ship pump for repair

## SPECIFICATIONS

### Vacuum pumps

<i>type</i>		ROTAVAC valve control	ROTAVAC vario control	ROTAVAC valve tec	ROTAVAC vario tec
max. suction power ( 50 / 60 Hz)	m <sup>3</sup> /h	1,7	1,6	0,75 / 0,9	1,0
final pressure (absolute)	mbar	9	2	12	12
final pressure (absolute) with gas-ballast	mbar	20	4	20	20
max. adm. counterpressure at output end (absolute)	bar	1,1			
adm. ambient temperature storage	°C	-10 to +60			
operation	°C	+10 to +40			
admissible relative humidity for operation (non-dewing)	%	30 to 85			
dissipated power	W	180	160	80	160
idle speed	1/min	1500	400 - 1500	1500 / 1800	400 - 1500
rated speed					
motor protection		Thermal winding protection			
sound pressure level (1 m distance)	dBA	42	<45	42	<45
hose shaft in / out	mm	8			
size: l x w x h	mm	235 x 145 x 175	236 x 143 x 205	315 x 145 x 190	236 x 143 x 180
weight, approx.	kg	10,5	6,5	6	5,9

### material in touch with media

components	ROTAVAC valve control	ROTAVAC vario control ROTAVAC vario tec ROTAVAC valve tec
cover, inside	PTFE, carbon fiber reinf.	PTFE
lid	ETFE	PFA
cocking disc, diaphragm	ETFE	PFA
valve	FFKM	FFKM
diaphragm	PTFE-NBR sandwich	PTFE
hoses	PTFE	PTFE
fittings / in / out	ETFE	ETFE

## VAC control automatic

feature	value
range	1400 - 1 mbar
resolution	1 mbar
accuracy	0,5% f.s.
materials in touch with media	ceramics, PP, Viton
measuring cycles	20 ms
operating temperature	0 ... 40°C
storage temperature	-10 ... +50°C
power supply	90 - 260 VAC 50/60 Hz
dissipated power	max. 13 W
switching outputs	24 VDC, max. 5W
terminal for vacuum	hose shaft, 9 mm
terminal for aeration	hose shaft 4mm
state of aeration valve powerless	open
display	LED 13mm
size	145x85x50(85) mm
weight	550g

## WARRANTY, LIABILITY & COPYRIGHT

### Warranty

Heidolph Instruments provides a three-year warranty on the products described here (with the exception of glass and consumable parts) if registered with enclosed warranty card or via internet ([www.heidolph.com](http://www.heidolph.com)). Warranty starts with the date of registration. Without registration warranty starts according to serial number. This warranty covers defects in materials and workmanship. Transit damage is excluded from this warranty.

To file for such warranty service, contact Heidolph Instruments (phone ++49-9122-9920-68) or your local Heidolph Instruments Dealer. If defects in material or workmanship are found, your item will be repaired or replaced at no charge.

Misuse, abuse, neglect or improper installation are not covered by this warranty.

Alterations to the present warranty need Heidolph Instruments' consent in writing.

### Exclusion Clause

Heidolph Instruments cannot be held liable for damage from improper use or misuse. Remedy for consequential damage is excluded.

### Copyright

Copyright in pictures and wording of the present Instruction Manual is held by Heidolph Instruments.

## FAQ / REPAIR WORK

If any **aspect** of installation, operation or maintenance remains unanswered in the present Manual, please contact the following address:

For repair services please call Heidolph Instruments (phone: +49 - 9122 - 9920-68) or your local, authorized Heidolph Instruments Dealer.

### Note:

You will receive approval for sending your defective item to the following address:

**Heidolph Instruments GmbH & Co. KG**  
**Lab Equipment Sales**  
**Walpersdorfer Str. 12**  
**D-91126 Schwabach / Germany**  
**phone: ++49-9122-9920-68**  
**Fax: ++49-9122-9920-65**  
**E-Mail: [Sales@Heidolph.de](mailto:Sales@Heidolph.de)**

### Safety Information

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact *substance*
- take proper *protective* measure to ensure the safety of our receiving and service personnel
- *mark* the pack IAW Hazardous Materials Act

## CE-DECLARATION OF CONFORMITY

We herewith declare that the present product complies with the following standards and harmonized documents:

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