
**METTLER TOLEDO**

**InPro 6850 i, 6900 i & 6950 i  
Series O<sub>2</sub> Sensors**

**Instruction manual  
Bedienungsanleitung  
Instructions d'utilisation**

**English** **Page 3**

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InDip, InFit, InPro, ISM, and InTrac are registered trademarks of the Mettler Toledo Group in Switzerland and a further twelve countries.

# **InPro 6850 i, 6900 i & 6950 i Series O<sub>2</sub> Sensors**

## **Instruction manual**

# Contents

<b>1</b>	<b>Introduction</b>	<b>5</b>
<b>2</b>	<b>Important notes</b>	<b>6</b>
2.1	Notes on operating instructions	6
2.2	Intended use	6
2.3	Safety instructions	7
2.4	Examples of some typical applications	8
2.5	Use in Ex-zones	8
2.6	Ex-classification ATEX	9
2.6.1	Introduction	9
2.6.2	Rated data	9
2.6.3	Special conditions	10
2.7	 Ex-classification FM approved	11
<b>3</b>	<b>Product description</b>	<b>12</b>
3.1	General information	12
3.2	Principle	12
3.3	Scope of delivery	13
3.4	Equipment features	14
<b>4</b>	<b>Installation</b>	<b>15</b>
4.1	Mounting the sensor	15
4.2	Connection	15
4.2.1	Connecting the sensor to a AK9 cable	15
4.2.2	Connecting the AK9 cable to the transmitter	16
<b>5</b>	<b>Operation</b>	<b>17</b>
5.1	Start-up and polarizing	17
5.2	Calibration	18
5.2.1	Purpose of calibration	18
5.2.2	What you have to know for calibration	18
5.2.3	Single point calibration	19
5.2.4	Dual point calibration	19
<b>6</b>	<b>Maintenance</b>	<b>20</b>
6.1	Inspection of the sensor	20
6.1.1	Visual inspection	20
6.1.2	Testing the METTLER TOLEDO O <sub>2</sub> sensor master	21
6.1.3	Testing the sensor via a transmitter	22
6.1.4	ISM design	23
6.2	Changing the electrolyte, the membrane body or the interior body	23
<b>7</b>	<b>Storage</b>	<b>26</b>
<b>8</b>	<b>Product specification</b>	<b>26</b>
8.1	Certificates	26
8.2	Specifications	27
<b>9</b>	<b>Ordering information</b>	<b>29</b>
9.1	Sensors with ISM functionality	29
9.2	Accessories	29
9.3	Spare parts	29
9.4	Recommended transmitters	30
9.5	Recommended housings	31
<b>10</b>	<b>Theory of the polarographic sensor</b>	<b>32</b>
10.1	Introduction	32
10.2	Principle of the design of an oxygen electrode	32
10.3	Parameters determining current	36
10.4	Polarization voltage	36
10.5	Temperature	37
10.6	Dependence on flow	37
10.7	Oxygen partial pressure – oxygen concentration	37

## 1 Introduction

Thank you for buying the **InPro® 6850 i/6900 i/6950 i sensor from METTLER TOLEDO.**

The construction of the InPro series employs leading edge technology and complies with safety regulations currently in force. Notwithstanding this, improper use could lead to hazards for the user or a third-party, and/or adverse effects on the plant or other equipment.

**Therefore, the operating instructions must be read and understood by the persons involved before work is started with the sensor.**

The instruction manual must always be stored close at hand, in a place accessible to all people working with the InPro sensor.

If you have questions, which are not or insufficiently answered in this instruction manual, please contact your METTLER TOLEDO supplier. He will be glad to assist you.

## 2 Important notes

### 2.1 Notes on operating instructions

These operating instructions contain all the information needed for safe and proper use of the InPro 6850i / 6900i / 6950i sensor.

The operating instructions are intended for personnel entrusted with the operation and maintenance of the sensors. It is assumed that these persons are familiar with the equipment in which the sensor is installed.

#### Warning notices and symbols

This instruction manual identifies safety instructions and additional information by means of the following symbols:

**This symbol draws attention to safety instructions and warnings of potential danger which, if neglected, could result in injury to persons and/or damage to property.**

**This symbol identifies additional information and instructions which, if neglected, could lead to defects, inefficient operation and possible loss of production.**

### 2.2 Intended use

**METTLER TOLEDO InPro 6850i/6900i/6950i sensors are intended solely for inline measurement of the oxygen partial pressure in liquids and gases, as described in this instruction manual.**

Any use of these sensors which differs from or exceeds the scope of use described in this instruction manual will be regarded as inappropriate and incompatible with the intended purpose.

The manufacturer/supplier accepts no responsibility whatsoever for any damage resulting from such improper use. The risk is borne entirely by the user/operator.

Other prerequisites for appropriate use include:

- compliance with the instructions, notes and requirements set out in this instruction manual.
- acceptance of responsibility for regular inspection, maintenance and functional testing of all associated components, also including compliance with local operational and plant safety regulations.
- compliance with all information and warnings given in the documentation relating to the products used in conjunction with the sensor ( housings, transmitters, etc).
- observance of all safety regulations governing the equipment in which the sensor is installed.

- correct equipment operation in conformance with the prescribed environmental and operational conditions, and admissible installation positions.
- consultation with Mettler-Toledo Process Analytics in the event of any uncertainties.

## 2.3 Safety instructions

- The plant operator must be fully aware of the potential risks and hazards attached to operation of the particular process or plant. The operator is responsible for correct training of the workforce, for signs and markings indicating sources of possible danger, and for the selection of appropriate, state-of-the-art instrumentation.
- It is essential that personnel involved in the commissioning, operation or maintenance of these sensors or of any of the associated equipment (e.g. housings, transmitters, etc.) be properly trained in the process itself, as well as in the use and handling of the associated equipment. This includes having read and understood this instruction manual.
- The safety of personnel as well as of the plant itself is ultimately the responsibility of the plant operator. This applies in particular in the case of plants operating in hazardous zones.
- The oxygen sensors and associated components have no effect on the process itself and cannot influence it in the sense of any form of control system.
- Maintenance and service intervals and schedules depend on the application conditions, composition of the sample media, plant equipment and significance of the safety control features of the measuring system. Processes vary considerably, so that schedules, where such are specified, can only be regarded as tentative and must in any case be individually established and verified by the plant operator.
- Where specific safeguards such as locks, labels, or redundant measuring systems are necessary, these must be provided by the plant operator.
- A defective sensor must neither be installed nor put into service.
- Only maintenance work described in this operating instruction may be performed on the sensors.
- When changing faulty components, use only original spare parts obtainable from your METTLER TOLEDO supplier (see spare parts list, "Section 9.3").
- No modifications to the sensors and the accessories are allowed. The manufacturer accepts no responsibility for damages caused by unauthorised modifications. The risk is borne entirely by the user.

## 2.4 Examples of some typical applications

Below is a list of examples of typical fields of application for the oxygen sensors. This list is not exhaustive.

### Measurement in liquids:

- Biotech
- Chemical applications
- Brewing
- Beverage filtration
- Filling stations

### Measurement in gases:

- CO<sub>2</sub> recovery
- CO<sub>2</sub> purity
- Product storage
- Inert production

## 2.5 Use in Ex-zones

**Attention!** For an installation in Ex-zones please read the guidelines following hereafter:

 1258

### Ex-classification ATEX:

Ex ia IIC T6/T5/T4/T3 Ga/Gb

Ex ia IIIC T69°C/T81°C/T109°C/T161°C Da/Db

### Number of the test certificate:

SEV 14 ATEX 0169 X  
IECEX SEV 14.0026X

### Ex-classification FM approved:

IS/I, II, III/1/ABCDEFGH/T6 Ta = 60°C  
– 53 800 002; Entity

## 2.6 Ex-classification ATEX

### 2.6.1 Introduction

According to Directive 94/9/EC (ATEX 95) Appendix I, the O<sub>2</sub> oxygen sensors type InPro 6XXX is a devices of equipment group II, category 1/2G which, according to Directive 99/92/EC (ATEX 137) can be used in zones 0/1 or 1/2 or 1 or 2 as well as gas groups IIA, IIB and IIC, which are potentially explosive due to combustible substances in the temperature T3 to T6.

The requirements specified in EN 60079-14 must be observed during use/installation.

According to Directive 94/9/EC (ATEX 95) Appendix I, the O<sub>2</sub> oxygen sensors type InPro 6XXX is a devices of equipment group II, category 1/2D which, according to Directive 99/92/EC (ATEX 137) can be used in zones 20/21 or 21/22 or 21 or 22, which are potentially explosive due to combustible dust.

The requirements specified in EN 60079-14 must be observed during use/installation.

For the analog version of the O<sub>2</sub> oxygen electrode, the O<sub>2</sub> measurement circuit, temperature measurement circuit and data chip circuit are part of the common intrinsically safe system and are jointly connected to and operated by a separately certified transmitter.

The digital version of the O<sub>2</sub> oxygen sensor is connected to and operated by two-wire cable to the certified transmitter.

The intrinsically safe circuits are galvanically isolated from the non-intrinsically safe circuits up to a nominal voltage peak value of 375 V and from the earthed parts up to a nominal voltage peak value of 30 V.

### 2.6.2 Rated data

#### **Analog O<sub>2</sub> Oxygen sensor**

With type of protection: intrinsic safety to Ex ia IIC

#### **O<sub>2</sub> measuring circuit, temperature measuring circuit and data chip circuit**

Only for connection to certified intrinsically safe circuits. Maximum values:

$$U_i \leq 16 \text{ V}, I_i \leq 190 \text{ mA}, P_i \leq 200 \text{ mW}$$

$$L_i = 0 \text{ (effective internal inductance)}$$

$$C_i = 900 \text{ pF (effective internal capacitance)}$$

The values above apply, each as the sum of all the individual circuits of the associated intrinsically safe supply and evaluation unit (transmitter).

#### **Digital O<sub>2</sub> Oxygen sensor**

With type of protection: intrinsic safety to Ex ia IIC

**Two-wire current circuit**

Only for connection to certified intrinsically safe circuits. Maximum values:

$$U_i \leq 16 \text{ V}, I_i \leq 30 \text{ mA}, P_i \leq 50 \text{ mW}$$

$$L_i = \text{negligible}$$

$$C_i = \text{negligible}$$

**2.6.3 Special conditions**

- The relationship between the maximum permissible ambient or media temperature and temperature class, for category 1G applications, zone 0, is shown in the following table:

<b>Temperature class</b>	<b>Max. ambient or media temperature</b>
T 6	68 °C
T 5	80 °C
T 4	108 °C
T 3	160 °C

- The relationship between the maximum permissible ambient or media temperature and temperature class, for category 1D applications, zone 20, is shown in the following table:


<b>Temperature class</b>	<b>Max. ambient or media temperature</b>
T 69 °C	68 °C
T 81 °C	80 °C
T 109 °C	108 °C
T 161 °C	160 °C

- The capacitance and inductance of the connecting cable has to be considered.
- The O<sub>2</sub> Oxygen sensor type InPro 6XXX can be used in/with the fittings InFit<sup>®</sup> 76\*-\*<sup>\*\*\*</sup> or InTrac<sup>®</sup> 7\*\*-\*<sup>\*\*\*</sup>, or in/with other suitable fittings in potentially explosive areas.
- The metal body of the O<sub>2</sub> Oxygen sensors, or the fittings InFit 76\*-\*<sup>\*\*\*</sup> or InTrac 7\*\*-\*<sup>\*\*\*</sup>, or other appropriate fitting is optionally included in the routine pressure test of the system.
- The independent fitting used for installation of O<sub>2</sub> Oxygen sensor must be conductively connected to the equipotential bonding system.

2.7

Ex-classification FM approved

**HAZARDOUS (CLASSIFIED) LOCATION**

Class I, Division 1, Groups A, B, C and D  
 Class II, Division 1, Groups E, F, and G  
 Class III, Division 1  
 T6 Ta=60 °C

Entity Parameters:  
 Vmax (Uij)=16V, Imax (Ii)=50mA, Pi=0.25W  
 Ci=0.1 µF, Li=0 mH

**NONHAZARDOUS LOCATION**

Any FMRC Approved Single Multi-Channel Barrier or Apparatus

- Notes:
1. No revision to this drawing is permitted without FMRC approval.
  2.  $V_{max} > V_t$ ;  $I_{max} > I_t$ ;  $(C_i \text{ of all loops} + C \text{ cable}) < C_a$ ;  $(L_i \text{ of all loops} + L \text{ cable}) < L_a$   
 $P_{max} \text{ or } P_i > P_o$
  3. Single Multi-Channel IS Barrier or Apparatus must be FMRC Approved.
  4. Single Multi-Channel IS Barrier or Apparatus manufacturer's control drawings must be followed when installing the System. IS Barrier or Equipment may be installed within the Hazardous (Classified) location for which it is approved.
  5. Installation must be in accordance with Article 500 of the NEC®(ANSI/NFPA 70) and ANSI/ISA RP-12.6.

WARNING: substitution of components may impair intrinsic safety.

<b>METTLER TOLEDO</b>	Process Analytics CH-8802 Udorf	Drawn F. Kogelmann 11-Dez-03	Scale 1:1
Schutzmerk ISO 16016 beachten / Refer to protection notice ISO 16016		Changed Dario Meier 25-Mai-16	Format DIN A3
Title <b>FM control drawing</b>		Note -	Replaces -
Model: MODEL NAME E.4	ArtikelNr: 53800002	Drawing No <b>53800002</b>	Index <b>E.4</b>
		Drawing: E.4	RelLevel: Released
			Sheet: 1/1

## 3 Product description

### 3.1 General information

The **oxygen sensor series InPro 6850 i/6900 i/6950 i** with integrated temperature probe are used for measurement of oxygen at low and medium concentrations.

The sensors are **sterilizable**, and in the majority of cases autoclavable too and they are **compatible with CIP** (Cleaning-In-Place).

InPro 6xxx i sensors with ISM<sup>®</sup> functionality offer Plug and Measure as well as enhanced diagnostics features.

### 3.2 Principle

#### **Amperometric oxygen sensors:**

The amperometric oxygen sensors of METTLER TOLEDO base all on the same measuring principle according to Clark. However different series are offered, which clearly differ in the number and arrangement of their electrodes and thus in their specifications:

- The measurement system of the InPro sensors consists of a working electrode (cathode), a counter electrode (anode), a reference electrode and a guard ring electrode. The measurement system is separated from the process medium by an oxygen permeable membrane.
- The measuring technique in principle is alike with all sensors.
- The transmitter supplies a constant voltage between cathode and anode.
- The electrolyte creates a conductive connection between the electrodes.
- The oxygen molecules migrate from the measurement medium through the membrane to the cathode to which the voltage is supplied and are reduced. At the same time, oxidation takes place at the anode.
- This causes a current to flow between the anode and cathode which is directly proportional to the partial pressure of oxygen (pO<sub>2</sub>) in the process medium. With digital ISM sensors, the sensor itself converts the current into oxygen concentration and communicates the value to the transmitter.
- The guard ring electrode of the InPro 6900 i and 6950 i built into the sensor reduces the oxygen that migrates to the side of the cathode and can distort the measurement. The guard ring electrode therefore enables traces of oxygen to be precisely determined at even the lowest concentrations.

**Note: Please refer to “Section 10 – Theory of the polarographic sensor” for further information.**

**ISM sensors:**

All oxygen sensors with the index "i" (6850 i, 6900 i, 6950 i) are equipped with ISM.

Principle: In the sensor head a chip is integrated, which takes over the entire monitoring and control of the sensor and, beyond that, stores all sensor data. This chip is responded via the transmitter.

The following data are available permanently in the sensor:

- type of sensor
- serial no.
- software version
- hardware version
- order no.
- operating time
- calibration time and calibration date
- calibration table

To check the system, the following indicators are supervised:

- temperature
- slope
- zero current
- air current
- polarization voltage

Based on these informations, the transmitter calculates a wear monitor and displays it depending upon transmitter type differently. (see the respective instruction manuals)

ISM enables the connection of the digital sensor to iSense. This user-friendly software allows to manage all sensor information and store it in a data base. Furthermore, the sensor can be calibrated via the iSense software.

### 3.3 Scope of delivery

Each sensor is supplied fully assembled and factory-tested for correct function together with:


- a quality control certificate
- inspection certificates 3.1  
(complying with EN 10204)

METTLER TOLEDO ISM O<sub>2</sub> sensors are supplied with fitted membrane body but without electrolyte, covered with the protection cap. They have been checked for proper function. See electrolyte order details at section „Spare Parts“, 9.3.

Digital sensors must be filled with electrolyte before start-up.

### 3.4 Equipment features

#### 12 mm sensor

METTLER TOLEDO ISM O<sub>2</sub> sensors are supplied with fitted membrane body but without electrolyte and covered with the protection cap and have been checked for proper function.

## 4 Installation

### 4.1 Mounting the sensor


**Important! Remove the protection cap before mounting the sensor.**

#### **Mounting the sensor in a housing**

Please refer to the instruction manual of your housing explaining on how to mount the sensor in place.

#### **Mounting the sensor directly on a pipe or a vessel**

The 12 mm sensors can be mounted directly through a socket with inside thread Pg 13.5 and securely tightened via the Pg 13.5 threaded sleeve.

### 4.2 Connection

#### 4.2.1 Connecting the sensor to a AK9 cable



The sensor is connected to the transmitter via a AK9 cable. The AK9 cable ensures a secure connection between the transmitter and the sensor under harsh industrial conditions. The robust watertight IP 68 connector housing guarantees maximum process safety.


**Do not touch the sensor at the AK9 connector plug!**

Tightly screw the plug to fasten the two parts.

#### 4.2.2 Connecting the AK9 cable to the transmitter

**Note:** For connecting the AK9 cable to the terminals of the transmitter, please refer to the instructions given in the METTLER TOLEDO transmitter manual.

## 5 Operation

**Important! Before using the sensors for the first time, the electrolyte must be filled in (see "Chapter 6.2").**

### 5.1 Start-up and polarizing

**Important! The protection cap must be removed before mounting the sensor in the process.**

When the system is operated for the first time or if the sensor has been disconnected from the voltage source (transmitter or O<sub>2</sub> sensor master) for longer than 5 minutes, the sensor has to be polarized prior to calibration by connecting it to the operating O<sub>2</sub> transmitter or to a sensor master. After 6 hours, the sensor is fully polarized and ready for operation.

**Attention: The InPro 6950 i should never be polarized at air!**

During this time, it is recommended to keep the protection cap filled with fresh conditioning solution on the sensor (see "Section 7 – Storage"), especially if the sensor is being polarized for longer than 6 hours.

A shorter polarization period is sufficient if the sensor has been disconnected for only a few minutes. The following table serves to establish the correct polarization time in relation to the depolarization time.

<b>Depolarization time<sup>1</sup> t<sub>depol</sub> [Min.]</b>	<b>Minimum required polarization time<sup>2</sup> [Min.]</b>
> 30	360
15 to 30	6 × t <sub>depol</sub>
5 to 15	4 × t <sub>depol</sub>
< 5	2 × t <sub>depol</sub>

- 1** Depolarization time: Time span in which the polarization voltage is cut off from the sensor. This is the case:
  - during the time the cable is disconnected or no transmitter or sensor master is connected to the cable, or the transmitter has been disconnected from the current supply.
  - after changing the electrolyte or membrane body. In this case at least 6 hrs of polarization must follow.
- 2** Polarization time: Time span during which the sensor is under a polarization voltage.

**Important! Setting of the polarization voltage on the transmitter for correct measurements:**

- Standard applications InPro 6850 i: –675 mV

- Measurement of permanently low oxygen concentrations (< 500 ppb in liquids or < 10,000 ppm [vol.] in gases) in the presence of volatile acidic components (e.g. carbon dioxide during measurements in breweries) e.g. InPro 6900i/6950i: **-500 mV**

**Note: To ensure the supply of the correct polarization voltage the transmitter must be set accordingly.**

## 5.2 Calibration

### 5.2.1 Purpose of calibration

Each oxygen sensor has its own individual slope and own individual zero point. Both values are subject to change, for example, through electrolyte consumption or after exchange of electrolyte or membrane body. To ensure high measurement accuracy of the sensor, a calibration must be carried out regularly, but at least after each change of electrolyte or membrane. Prior to calibration, the sensor has to be polarized for at least 6 hours.

**Please remove the protection cap from the sensor, rinse the sensor with water and dry it for at least 10 minutes before starting calibration.**

**To check if your sensor needs a recalibration, you may dry it and take it in the air to check that the reading is close to 100 %. If not, then the sensor needs a new calibration.**

With ISM sensors all calibration data are stored in the sensor.

### 5.2.2 What you have to know for calibration

Calibration is generally carried out at the polarization voltage at which measurements are taken.

Zero-point calibration is recommended:

- when measuring less than 5 ppb in liquids
- when measuring less than 125 ppm (by vol.) in gases

**General remarks:**

- **For calibration in air, the sensor membrane must be dry**, since adhering water drops can falsify the measured oxygen value.
- Make sure that the **oxygen saturation index** of the calibration medium is **correct and remains constant** during calibration.
- In the event of calibration in water or sample medium, the **calibration medium must be in equilibrium with the air**. Oxygen exchange between water and air is only very slow. Therefore it takes quite long time until water is saturated with atmospheric oxygen.

- For correct calibration, a minimum flow rate of the calibration medium is necessary.
- Make sure that all other parameters, such as temperature and pressure, are constant.

For continuous applications, we recommend **periodic recalibration in line with your requirements on accuracy, the type of process in operation and your own experience**. The frequency of the need for re-calibration depends very much on the specific application, and therefore appropriate intervals cannot be exactly defined here.

### 5.2.3 Single point calibration

By carrying out a single point calibration, the factual slope of the sensor can be established. The calibration medium can be water with known oxygen saturation index (e.g. air-saturated water) or air with known water-vapor saturation (e.g. water-vapor saturated air).

After the sensor signal has stabilized, the complete measuring system can then be calibrated to the 100 % value of the desired measurable variable, e.g. 100 % air, 20.95 % O<sub>2</sub>, or 8.26 ppm at 25 °C (77 °F) and normal pressure (see instruction manual for the transmitter).

### 5.2.4 Dual point calibration

By carrying out a dual point calibration both slope and zero point of the sensor can be established.

**Important! In case of a dual point calibration, always start by the zero point calibration before calibrating the slope.**

As a rule, the zero point should be adjusted to zero manually, or it is automatically performed by the transmitter (see instruction manual for the transmitter).

Zero-point calibration is recommended:

- when measuring less than 5 ppb in liquids
- when measuring less than 125 ppm (by vol.) in gases

**Attention! Incorrect zero point calibration is a frequent source of measurement error. For correct calibration, we recommend the use of carbon dioxide gas with a level of purity of at least 99.9995 %.**

After the sensor signal has stabilized (after **6 ... 12 hours**), the sensor can be calibrated through the relevant transmitter to the zero % value of the desired measurable variable, e.g. 0 % air, 0.0 % O<sub>2</sub>, or 0.0 ppm (see instruction manual for the transmitter).

## 6 Maintenance

### 6.1 Inspection of the sensor

#### 6.1.1 Visual inspection

To check your sensor, we recommend the following procedure:

- The contacts of the connector must be dry. Moisture, corrosion and dirt in the connector can lead to false readings.
- Check the cable for buckling, brittle areas or ruptures.
- Before calibration always examine the membrane foil optically for signs of damage. The foil must be intact and clean. Dirty membranes should be wiped clean using a soft, moist tissue.

**Note: An undulated membrane has no influence on the sensor performance, assuming the membrane is intact.**

- The membrane body must be replaced if the sensor has too long a response time, the reading is unstable or subject to drift, and if the sensor cannot be calibrated or the membrane shows sign of mechanical damage.
- Check the cathode area for discoloration, contamination or cracks in the glass. If necessary rinse with demineralized water and clean with a clean soft brush or soft paper tissue.

**Attention! Do not use any cleaning agents containing alcohol. This could damage the sensor or lead to fault current.**


**Attention! The glass body is fragile and sensitive to vibration.**

### 6.1.2 Testing the METTLER TOLEDO O<sub>2</sub> sensor master

We recommend to use the METTLER TOLEDO O<sub>2</sub> sensor master to check the quality of your sensor as follows:

- **Connect the sensor to the O<sub>2</sub> sensor master.**

As soon as the sensor is connected to the O<sub>2</sub> sensor master, the polarization function is automatically activated. Please note: if the sensor was disconnected from the transmitter for longer than 5 minutes, the sensor must be polarized first (polarizing time see "Section 5.1") to get representative test results.

- **Battery Check:**  
Push the switch to the left. If the battery is ok and the O<sub>2</sub> sensor master is operational the green LED lights up. Otherwise, please consult the instruction manual of the O<sub>2</sub> sensor master.
- **Sensor Check:**  
For this test the O<sub>2</sub> sensor must be fully polarized and the membrane of the sensor must be dry and clean.

Expose the sensor (connected to O<sub>2</sub> sensor master) to the air. By pushing the switch to the right to the position "2 – Sensor check", the O<sub>2</sub> sensor master checks whether the electrode current for air measurement delivered by the sensor is within the admissible range, i.e. 2500 to 6000 nA for InPro 6950.

If the green LED lights up the current for air measurement is within the admissible range.

If the LED does not light up, you should check the battery of the O<sub>2</sub> sensor master (see instruction manual "Accessories"). If the battery is working, than there is probably a problem with your sensor. You should change the electrolyte and/or the membrane body of your sensor. If after a membrane change the LED still does not light up, this means that there is maybe something wrong with the interior body of the sensor. You should then change it (see "Section 6.2").

**Important! The Sensor Check function only verifies the correctness of the electrode current for air measurement. In order to be absolutely sure of the functionality of the sensor, the residual signal in an oxygen free medium should also be controlled (see "Section 6.1.3").**

### 6.1.3 Testing the sensor via a transmitter

A periodic zero current measurement (**no zero point calibration!**) is recommended for verification of proper sensor function.

**Note: At the time you carry out the zero current measurement, the sensor must be polarized.**

Zero current measurement can be done by using zeroing gel (order no. 30 300 435) or nitrogen (N<sub>2</sub>) or carbon dioxide (CO<sub>2</sub>) calibration gases with a purity of at least 99.995 %, alternatively in a sample medium saturated with one of these gases.

After 2 minutes in an oxygen-free sample medium, the reading on the transmitter should drop to below 10 % of the reading in ambient air, and within 10 minutes the value should have dropped to below 1 %.

If the measured values are too high, this suggests a depleted electrolyte or a defective membrane. In the first instance replace the electrolyte, and in the second case exchange both the membrane body and the electrolyte accordingly.

If after such procedures the above mentioned values are still not reached, replace the interior body. If this doesn't solve the problem too send the sensor to your local METTLER TOLEDO representative for inspection.

Many sample media contain volatile substances which, even at very low concentrations, have a clearly perceptible smell. Similarly to oxygen, these substances are able to invade the electrolyte through the gas-permeable membrane. Accordingly, they become noticeable when changing the electrolyte. In most cases, such substances have absolutely no

influence on the measuring properties of the sensor. Slight discoloration of the electrolyte also has no effect on the measuring properties.

#### 6.1.4 ISM design

The integrated ISM functionality allows an extensive monitoring of the sensor. The following parameters are stored in the sensor:

- serial no.
- type of sensor
- order no.
- calibration data
- CIP/SIP counter
- slope
- zero point

When starting the following automatic test procedures are implemented:

- digital communication
- plug & measure
- pre-calibration
- predictive maintenance

## 6.2 Changing the electrolyte, the membrane body or the interior body

**Note: The InPro 6900 i and 6950 i use a special electrolyte which contains an oxygen scavenger. This electrolyte ensures fast response time and together with the guard ring increases the signal stability of the sensor. The electrolyte must be exchanged on a regular basis or when the sensor has been exposed to the air for more than 24 hours without a watering cap filled with conditioning solution.**

**Note: The content of the electrolyte bottle has to be used within 3 months after opening.**

If the membrane and/or the interior body exhibits signs of failure (long response time, increased current in an oxygen-free medium, mechanical damage, etc.) the membrane body and/or the interior body has to be replaced.

**Warning! The O<sub>2</sub> electrolyte has an alkaline pH value. Contact of electrolyte with mucous membrane or eyes is to be avoided. Therefore protective gloves and safety glasses have to be worn for the following dismantling works. If such contact occurs, the affected area should be well rinsed with water. In the case of accident, or should ever any adverse signs appear, get immediate medical attention.**

When changing the electrolyte, the membrane body or the interior body, please observe the following instructions (see also the following illustration):

**Attention! Make sure that this maintenance step is carried out in clean place.**

1. Unscrew the cap sleeve from the sensor shaft and carefully pull it off the sensor.
2. Pull off the membrane body from the interior body. If it is tight-fitted, eject by pushing it with the flat finger tip. Before electrolyte is refilled, the membrane body must be removed from the cap sleeve!
3. Rinse the interior body with demineralized water and carefully dab it dry with a paper tissue.

**Note: steps 4 to 7 may only be carried out when changing the interior body.**

4. Unscrew the retainer nut of the interior body with an adjustable wrench or with a 9 mm wrench.
5. Remove the interior body by pulling it out of the sensor shaft. If necessary use a plier.


**Warning! Do not twist the interior body. Otherwise the connection pins can be damaged.**

6. Insert the new interior body in the sensor shaft. Turn the interior body in the shaft until the slit of the interior body is aligned with the pin placed in the shaft.
7. Press the body in the shaft and screw the new retainer nut in place.
8. Examine the O-rings visually for mechanical defects, and replace if necessary.
9. Half-fill the new membrane body with O<sub>2</sub> electrolyte.

**Note: The electrolyte bottle is equipped with a special pouring system. To ensure proper functioning, hold the bottle vertically, upside-down.**


 **Note: make sure that all air bubbles are removed from the membrane body. Air bubbles can be removed by carefully tapping on the membrane body.**

10. Slip the membrane body over the interior body while holding the sensor in a vertical position. The excess electrolyte will be displaced and have to be removed with a paper tissue.

 **Important! No electrolyte, sample media or contamination may be present between the membrane body and the cap sleeve. Please check carefully!**

11. Carefully slip the cap sleeve over the fitted membrane body, holding the sensor in a vertical position and screw it tight. The cap sleeve must be clean and dry.
12. After each exchange of electrolyte or membrane body, the sensor has to be repolarized and recalibrated.

### Replacement of the interior body

## 7 Storage

For storage periods longer than 24 hours, we recommend you to use a protection cap filled with the same cleaning and conditioning solution as that used for our portable DO measurement system (order no. 52 200 255). This solution contains an oxygen scavenger which ensures that the electrolyte of the sensor does not come into contact with oxygen when not in use.

When the sensor is stored without polarization for more than one week, the electrolyte has to be removed.

Prepare the cleaning and conditioning solution as follows: Dissolve one tablet in 40 ml of deionized water and wait 5 minutes for the tablet to be completely dissolved. Fill the protection cap with this solution and place it over the tip of the sensor. This solution has some cleaning properties which will keep the membrane free of microorganisms. In case you do not have any cleaning and conditioning set, you may also use checking gel or deaerated water in the protection cap. Before mounting the sensor in the process, always remove the protection cap and rinse the tip of the sensor with water.

**Attention! If the storage period of the sensor without current supply (transmitter, sensor master) exceeds 1 week, the sensor should be stored dry, i.e. without any electrolyte in the membrane body. A sensor being stored dry (without electrolyte in the membrane body) may on no account be connected to the O<sub>2</sub> sensor master or any other polarization module.**

## 8 Product specification

### 8.1 Certificates

Each sensor is delivered with a set of **3.1 certificates** (complying with EN 10204).

All wetted metal parts (sensor shaft, cap sleeve and membrane body) are identified with an engraved symbol corresponding to the heat number on the paper certificate delivered with the sensor.

Each wetted metal part (sensor shaft, cap sleeve and membrane body) is polished in order to get a surface roughness lower than 0.4 µm (16 µin). This represents a roughness grade number of N5 (according to ISO 1320:1992).

## 8.2 Specifications

<b>InPro 6850i/6900i/6950i</b>	
<b>Measurement principle</b>	Polarographic Clark electrode
<b>Working conditions</b>	
Pressure resistance (measurement)	6850i: 0.2 ... 6 bar 6900i: 0.2 ... 6 bar (9 bar with T-6900R) 6950i: 0.2 ... 6 bar
Mechanical pressure resistance	Max. 12 bar
Temperature range (measurement)	0 ... 80 °C (32 ... 176 °F)
Temperature range (mechanical)	6850i: -5 ... 140 °C (23 ... 284 °F) 6900i: -5 ... 121 °C (23 ... 249.8 °F) [sterilizable and autoclavable] 6950i: -5 ... 121 °C (23 ... 249.8 °F) [sterilizable]
<b>Construction</b>	
Temperature compensation	Automatic
Cable connection	AK9
O-ring material	Silicone FDA and USP Class VI approved
Membrane material	PTFE/Silicone/PTFE (reinforced with steel mesh)
Material sensor body (wetted parts)	316L stainless steel
Surface roughness of wetted metal parts	N5 (R <sub>a</sub> = 0.4 µm [16 µin])
Quick disconnect interior body	Standard
Cathode	Pt
Anode	6850i: Pt 6900i: Ag 6950i: Pt
Guard ring	6850i: - 6900i: Pt 6950i: Pt
Reference	Ag
<b>Dimensions</b>	
Sensor diameter	6850i: 12/25 mm 6900i: 12 mm 6950i: 12 mm
Immersion length (a) for 12 mm sensor	6850i/6900i: 70, 120, 220, 320, 420 mm (2.8, 4.7, 8.66, 12.6, 16.5") 6950i: 70, 120, 220, 320 mm (2.8, 4.7, 8.66, 12.6")
Immersion length (a) for 25 mm sensor	6850i: 80, 160, 260, 360 mm (3.2, 6.3, 10.2, 14.2") 6900i: - 6950i: -
<b>Performance</b>	
Detection limit	6850i: 6 ppb to saturation 6900i: 1 ppb to saturation 6950i: 0.1 ppb to saturation
Accuracy	6850i: ≤ ± (1 % + 6 ppb) 6900i: ≤ ± (1 % + 1 ppb) 6950i: ≤ ± (1 % + 0.1 ppb)
Response time at 25 °C/77 °F	98 % of final value in < 90 s

Sensor signal in ambient air at 25 °C/77 °F	6850i: 50 ... 110 nA 6900i: 250 ... 500 nA 6950i: 2500 ... 5000 nA
Residual signal in oxygen-free medium	6850i: < 0.1 % of the signal 6900i: < 0.3 % of the signal 6950i: < 0.025 % of the signal

### Certification

EHEDG, 3A	Yes
3.1 B (EN 10204.3/1.B)	Yes
ATEX Certificate	6850i/6900i: Yes 6950i: Yes
FM Approval	6850i/6900i: Yes 6950i: Yes
FDA/USP VI	Yes
Quality Control	Yes

### Compatibility

with METTLER TOLEDO transmitters	see "Section 9.4"
with METTLER TOLEDO housings	see "Section 9.5"

## 9 Ordering information

For more detailed information refer to the technical data sheet. Ask your local distributor.

### 9.1 Sensors with ISM functionality

#### Ordering Information

– straight version	6850i	6900i	6950i
70 mm	52 206 118	52 206 316	52 206 127
120 mm	52 206 119	52 206 317	52 206 128
220 mm	52 206 120	52 206 318	52 206 129
320 mm	52 206 121	52 206 319	52 206 130
420 mm	52 206 122	–	–
– 25 mm version			
80 mm	52 206 123	–	–
160 mm	52 206 124	–	–
260 mm	52 206 125	–	–
360 mm	52 206 126	–	–

### 9.2 Accessories

Accessories	Order No.
O <sub>2</sub> sensor master digital ISM	52 206 329
Cables	
AK9 coax cable with K8S connector 1 m	59 902 167
AK9 coax cable with K8S connector 3 m	59 902 193
AK9 coax cable with K8S connector 5 m	59 902 213
AK9 coax cable with K8S connector 10 m	59 902 230
AK9 coax cable with K8S connector 20 m	52 300 204

### 9.3 Spare parts

Spare parts	Order No.
– InPro 6850 i	
Membrane body (single), T-96	52 200 071
Membrane kit T-96 (4 membranes, 1 O-ring set silicone, 25 ml of electrolyte, wetted parts SS 316L)	52 200 024
Membrane bodies (16 pcs), T-96	52 206 114
O <sub>2</sub> Electrolyte pack (3 × 25 ml)	30 298 424
Interior body (replacement), InPro 6850i	52 206 347
– InPro 6900 i	
Reinforced membrane body (single), InPro 6900 (T-6900R)	52 201 108
Reinforced membrane kit InPro 6900 (T-6900 R)	52 201 109
InPro 6900 electrolyte pack (3 × 5 ml)	30 298 425
Interior body (replacement), InPro 6900	52 200 943
– InPro 6950 i	
Membrane body (single)	52 206 105
Membrane kit (4 membranes, 1 O-ring set, 2 × 5 ml of electrolyte)	52 206 106
InPro 6950 electrolyte pack (3 × 5 ml)	30 298 426
Interior body (replacement), InPro 6950	52 206 112

## 9.4 Recommended transmitters

<b>Transmitters M700</b>	<b>Order No.</b>
Transmitter base, coated, M700C	52 121 171
Transmitter base, coated, Ex, VPW, M700XC/VPW	52 121 172
Transmitter base, coated, Ex, 24V, M700XC/24V	52 121 173
Transmitter base, stainless steel, M700S	52 121 174
Transmitter base, stainless steel, Ex, VPW, M700XS/VPW	52 121 175
Transmitter base, stainless steel, Ex, 24V, M700XS/24V	52 121 176

<b>Transmitter modules M700</b>	<b>Order No.</b>
ISM O <sub>2</sub> module 4700i	52 121 263
ISM O <sub>2</sub> module, Ex, 4700iX	52 121 263
ISM O <sub>2</sub> module 4700i ppb	52 121 265
ISM O <sub>2</sub> module, Ex, 4700iX ppb	52 121 266
ISM O <sub>2</sub> module 4700i traces	52 121 295
ISM O <sub>2</sub> module, Ex, 4700iX traces	52 121 294

<b>Transmitters M300D</b>	<b>Order No.</b>
M300i, ¼DIN, single-channel, multi-parameter	52 121 354
M300i, ½DIN, single-channel, multi-parameter	52 121 355
M300i, ¼DIN, dual-channel, multi-parameter	52 121 356
M300i, ½DIN, dual-channel, multi-parameter	52 121 357

<b>Transmitters M400</b>	<b>Order No.</b>
M400, type 1	52 121 348
M400, type 2	52 121 349
M400, type 3	52 121 350

For more information consult the "Parameter Fit Guide" below.

### Parameter Fit Guide for Transmitter M400

	<b>M400</b>	<b>M400</b>	<b>M400</b>
<b>Conventional sensors</b>	<b>Type 1</b>	<b>Type 2</b>	<b>Type 3</b>
pH/ORP	•	•	•
Cond 2-electrode	•	•	•
Cond 4-electrode	•	•	•
DO ppm	–	•	•
DO ppb	–	–	•
<b>ISM sensors</b>			
pH/ORP	•	•	•
Cond 2-electrode	•	•	•
Cond 4-electrode	•	•	•
DO ppm	–	•	•
DO ppm optical	–	–	•
DO ppb	–	–	•
DO traces	–	–	•

## 9.5 Recommended housings

### Housing (12 mm $\varnothing$ )

---

#### Static housing

---

InFit 761 CIP

---

#### Retractable housings

---

InTrac 777 e

---

InTrac 797 e

---

#### Immersion housing

---

InDip<sup>®</sup> 550

---

**Note: The housings are available in different versions. Please contact your local METTLER TOLEDO sales organization to get the correct ordering information.**

## 10 Theory of the polarographic sensor

### 10.1 Introduction

Two types of electrodes are employed in analytical work: **potentiometric** and **amperometric** electrodes.


- **Potentiometric electrodes** develop a voltage generated by the activity of a particular ion. Examples of such electrodes are glass electrodes (like pH electrodes) and most ion-selective electrodes. Their individual potentials cannot be determined. The measurable quantity is the difference of potential between the measuring electrode and an inert reference electrode. The potential of the reference electrode must be constant.

All potentiometric electrodes are subject to Nernst's law and for this reason electrodes and measuring instruments are in most cases interchangeable. An important requirement of potentiometric measurements is the virtually currentless determination of the electrode voltage. During measurement no chemical reaction occurs and the solution remains in equilibrium.

- In the case of **amperometric electrodes**, such as the **oxygen electrode**, activity measurement is based on a current measurement.

**The conventional oxygen electrode consists of a cathode and an anode conductively connected by an electrolyte.** A suitable polarization voltage between the anode and the cathode selectively reduces the oxygen at the cathode.

**Cathode reaction:**

These chemical reactions result in an electric current which is proportional to the oxygen partial pressure ( $p_{\text{O}_2}$ ). The oxygen electrode reduces constantly oxygen. Thereby, the concentration of dissolved oxygen is reduced. By diffusion this oxygen is replaced. The viscosity and flow rate of the solution are therefore important parameters. The electrode current of an oxygen electrode is determined not only by the oxygen partial pressure but by many other electrode parameters. The electrode currents of different electrode types may differ by several powers of ten. For this reason oxygen electrodes and amplifiers cannot be freely interchanged.

### 10.2 Principle of the design of an oxygen electrode

There are two main types of oxygen electrodes:

- Electrodes **without** membrane
- Electrodes **with** gas-permeable membrane (Clark Principle)

The membrane electrode according to Clark is today most widely used. As compared to the electrode without membrane it possesses the following advantages:

- Oxygen measurement in gases and solutions
- No mutual contamination of electrode and solution
- No or little dependency on flow


In the case of the Clark electrode, geometrical configuration is very important. In particular, the thickness of the electrolyte film between the cathode and the membrane must be within narrow tolerances so as to ensure good linearity and a low zero current (current in a nitrogen atmosphere).

METTLER TOLEDO oxygen sensors are available in different designs:

### Type A, 2-electrodes system, InPro 6800

InPro 6800 for medium and high oxygen concentrations. Cathode and anode/reference. Anode and reference are united in a silver/silver chloride electrode.


The following equilibration reaction takes place at the anode resp. reference:

### Type B, 3-electrodes system, InPro 6850 i


InPro 6850 i for medium and high oxygen concentrations. The reference corresponds to the conventional silver/silver chloride anode. The anode is a platinum electrode and separated from the reference.


The following reaction takes place at the anode:

The reference is a silver/silver chloride electrode. After polarization there is a stable equilibrium between electrode surface and electrolyte. No net-reaction takes place.


Reference reaction in equilibrium:



### Type C, 3-electrodes system, InPro 6900 (i)

InPro 6900 and InPro 6900i. For oxygen measurements in the lower ppb range. Here, the anode and reference are united in a silver/silver chloride electrode (as with type A). The sensors are equipped with an additional guard ring around the cathode. Like the cathode with the anode, this forms a closed electric circuit, which prevents that oxygen diffuses from the side to the cathode and falsifies the measurement result.


In all sensors the function of the cathode is identical.

Between cathode and reference a constant voltage of – 500 resp. – 675 mV (polarization voltage) is applied. Thereby a reduction of O<sub>2</sub> into OH<sup>-</sup> at the cathode takes place.


The total equation thus is:


### Type D, 4-electrodes system, InPro 6950 (i)

The InPro 6950(i) is intended for measuring permanent low oxygen level down to 0.1 ppb. In terms of measurement principle it is a combination of the InPro 6850 i and InPro 6900(i). It consists of 4 electrodes. The anode (platinum) and reference (silver/silver chloride) are separated in two electrodes. A guard ring is placed around the cathode. The cathode has the highest diameter from all amperometric oxygen sensors from METTLER TOLEDO


The electrodes show the following reactions:

Reference: No current

Due to the larger surface of the cathode a higher current appears during measurement. This allows for measurement of oxygen levels below 1 ppb.

As a consequence of the reaction at the cathode and the anode, a current flows between cathode and anode (4 electrons per oxygen molecule), which is proportional to the oxygen partial pressure at the cathode.

The amperage of the current flow depends on the surface of the cathode. Typical values are:

	<b>Air current</b>	<b>Zero current in % of the air current</b>
6800	50 – 110 nA	< 0.1
6850 (i)	50 – 110 nA	< 0.1
6900 (i)	250 – 500 nA	< 0.03
6950 i	2500 – 5000 nA	< 0.025

This current is measured and converted and indicated in the transmitter into an oxygen value. In ISM sensors this computation takes place in the sensor and the transmitters displays this value.

### 10.3 Parameters determining current

The quantity of oxygen diffused in and the magnitude of the electrode current are influenced by the following parameters:

- Oxygen partial pressure of the solution
- Membrane material and thickness
- Size of cathode
- Polarization voltage
- Temperature
- Flow conditions in the solution


Fick's law gives the mathematical relationship between these parameters:

$$I = k \times D \times a \times A \times \frac{pO_2}{x}$$

- I = Electrode current
- k = Constant
- D = Coefficient of diffusion of O<sub>2</sub> in the membrane
- a = Solubility of O<sub>2</sub> in the membrane material
- A = Cathode surface
- pO<sub>2</sub> = Oxygen partial pressure in the solution
- x = Thickness of gas-permeable membrane

### 10.4 Polarization voltage

The voltage between the anode and the cathode is so selected that the oxygen is fully (> A, see polarogram) reduced while other gases are unaffected (< D). The ideal voltage for the Pt/Ag/AgCl or Pt/Pt/Ag/AgCl system is between –500 and –750 mV.

The polarization voltage should remain as constant as possible. Besides a constant source of voltage the following prerequisites must be met: The electrical resistance of electrolyte film must not exceed a particular value so as to avoid a voltage drop.

The anode must have a large surface so as to prevent polarization of the anode by the electrode current.

## 10.5 Temperature

The temperature dependence of the current passing through an oxygen electrode when referred to a constant O<sub>2</sub> partial pressure is determined mainly by the properties of the gas-permeable membrane.

## 10.6 Dependence on flow

With most oxygen electrodes the electrode current is smaller in stagnant solutions than in agitated ones. In consequence of the oxygen consumption of the electrode, a reduction of oxygen takes place outside the membrane in close proximity to the cathode. The missing oxygen is replaced by diffusion. If the electrode current is strong, the solution cannot fully restore the reduced oxygen by diffusion. This results in an electrode current weaker than would correspond to conditions in the solution. In agitated solutions the oxygen is conveyed to the surface of the membrane not only by diffusion but additionally by the flow (convection). In that case no oxygen impoverishment occurs at the membrane surface.

A high degree of flow dependence occurs mainly with large cathodes, thin and highly permeable membranes, i.e. where electrode currents are large.

The problem of flow dependence is often solved by prescribing a minimum flow rate.

In METTLER TOLEDO InPro 6950 electrodes, the thin PTFE membrane determining the electrode current (i.e. the actual measuring signal) is separated from the sample solution by a relatively thick silicone membrane. This latter is highly permeable to oxygen molecules and thus acts as an oxygen reservoir. The diffusion of oxygen out of the sample solution into the silicone membrane is spread over a wide area. Since this results in less oxygen being extracted from the sample solution per unit area, the double PTFE/silicone membrane forms an effective buffer against disturbances due to hydrodynamic flow.

This membrane coupled with the guard ring and the special electrolyte ensures excellent signal stability even when the hydrodynamic flow stops (as on a beer filling line for example).

## 10.7 Oxygen partial pressure – oxygen concentration

The electrode current is dependent on the partial pressure of oxygen and the oxygen permeability of the membrane. Conversion of partial pressure into concentration of oxygen depends on the measurement medium (measurement in liquids or gases).

### Measurement in liquids

When measuring in liquids, the concentration of oxygen also depends on the solubility of the oxygen

in the measurement medium. However, since this is not measured by the sensor current, the concentration of oxygen must be calculated in the transmitter. To do this, Henry's law is applied which states that the concentration of oxygen is proportional to the partial pressure of oxygen ( $pO_2$ ).

$$CI = pO_2 \times \alpha$$

$\alpha$  = Solubility factor

If " $\alpha$ " is constant, the oxygen concentration can be determined by means of the electrode. This applies at constant temperature and with dilute aqueous solutions such as drinking-water.

The solubility factor is strongly influenced not only by the temperature but also by the composition of the solution:

<b>Medium, saturated with air</b>	<b>Solubility at 20 °C (68 °F) and 760 mm Hg</b>
Water	9.2 mg O <sub>2</sub> /l
4 mol/l KCl	2 mg O <sub>2</sub> /l
50 % Methanol-water	21.9 mg O <sub>2</sub> /l

Although the solubilities vary widely, the oxygen electrode gives the same reading in all three solutions.

Thus, determination of the oxygen concentration is only possible with constant and known solubility factors " $\alpha$ ".

Solubility may be determined by a Winkler titration or the method developed by Käppeli and Fiechter.

### Measurement in gases

The concentration of oxygen when measuring in gases is always given as a proportion by volume of the gas composition. Common units are % (by vol.) ppm (by vol.).

They can simply be converted from one unit of measurement to the other.

Example:

The percent by volume of the composition of air is generally known. For example, air contains 20.95 % oxygen. This corresponds to 209,500 ppm (by vol.). (Conversion: ppm value = 10,000 × value in %)

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## 2.7

## Classification EX – FM Approved

NONHAZARDOUS LOCATION		HAZARDOUS (CLASSIFIED) LOCATION	
<p>Any FMRC Approved Single Multi-Channel Barrier or Apparatus</p>		<p>Class I, Division 1, Groups A, B, C and D Class II, Division 1, Groups E, F, and G Class III, Division 1 T6 Ta=60 °C</p>	
<p>Notes:</p> <ol style="list-style-type: none"> <li>No revision to this drawing is permitted without FMRC approval.</li> <li><math>V_{max} &gt; V_t</math>; <math>I_{max} &gt; I_t</math>; <math>(C_i \text{ of all loops} + C \text{ cable}) &lt; C_a</math>; <math>(L_i \text{ of all loops} + L \text{ cable}) &lt; L_a</math> <math>P_{max}</math> or <math>P_i &gt; P_o</math></li> <li>Single Multi-Channel IS Barrier or Apparatus must be FMRC Approved.</li> <li>Single Multi-Channel IS Barrier or Apparatus manufacturer's control drawings must be followed when installing the System. IS Barrier or Equipment may be installed within the Hazardous (Classified) location for which it is approved.</li> <li>Installation must be in accordance with Article 500 of the NEC<sup>®</sup>(ANSI/NFPA 70) and ANSI/ISA RP12.6.</li> </ol>		<p>Entity Parameters: <math>V_{max} (U_i) = 16V</math>, <math>I_{max} (I_i) = 50mA</math>, <math>P_i = 0.25W</math> <math>C_i = 0.1 \mu F</math>, <math>L_i = 0 mH</math></p>	
<p>WARNING: substitution of components may impair intrinsic safety.</p>		<p>Probe</p>	
		<p>Drawn: F. Kogelmann 11-Dec-03 Scale: 1:1 Change: Darío Meier 25-Mar-16 Format: DIN A3 Note: - Replaces: - Drawing No: 53800002 Index: E.4</p>	
		<p>Process Analytics CH-8802 Urdorf Schutzvermerk ISO 16016 beachten! / Refer to protection notice ISO 16016 Title: FM control drawing Model: MODEL NAME E.4 ArticleNr.: 53800002 Drawing: E.4 ReLevel: Released Sheet: 1/1</p>	

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**France**      **Mettler-Toledo Analyse Industrielle S.A.S.**
**Germany**      **Mettler-Toledo GmbH**
**Switzerland**      **Mettler-Toledo (Schweiz) GmbH**