**METTLER TOLEDO**

## **InPro 7000-VP Series**

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## 2-electrode conductivity sensors

### Instruction manual

Subject to technical changes without prior notice.

© It is strictly forbidden to reprint this instruction manual or any parts thereof without the written permission of Mettler-Toledo GmbH, Process Analytics, Im Hackacker 15, 8902 Urdorf, Switzerland. No section or excerpt whatsoever may be reproduced or with the assistance of electronic systems be edited, duplicated or distributed, in particular in the form of photocopies, photographs, magnetic media or other recording methods. All rights reserved, especially the right of duplication and translation as well as in regard to patent and registration rights.

InPro is a registered trademark of the Mettler Toledo Group in Switzerland and a further twelve countries.

Mettler-Toledo GmbH, 8606 Greifensee, Switzerland

# Contents

|           |  |           |
|-----------|--|-----------|
| <b>1</b>  | <b>Product description</b>                       | <b>5</b>  |
| 1.1       | Introduction                                     | 5         |
| 1.2       | Equipment and scope of delivery                  | 5         |
| 1.3       | Technical data                                   | 7         |
| <b>2</b>  | <b>Safety</b>                                    | <b>8</b>  |
| 2.1       | Application compatibility                        | 8         |
| 2.2       | Proper utilization                               | 9         |
| 2.3       | Safety measures                                  | 9         |
| 2.4       | Use in Ex-zones                                  | 10        |
| <b>3</b>  | <b>Initial start-up</b>                          | <b>12</b> |
| 3.1       | Installation                                     | 12        |
| 3.2       | Electrical connections                           | 13        |
| <b>4.</b> | <b>Maintenance and troubleshooting</b>           | <b>14</b> |
| 4.1       | Conditions of warranty                           | 14        |
| 4.2       | Maintenance                                      | 14        |
| 4.3       | Troubleshooting                                  | 15        |
| <b>5</b>  | <b>Removal from operation, storage, disposal</b> | <b>16</b> |
| 5.1       | Removal from operation                           | 16        |
| 5.2       | Storage  | 16        |
| 5.3       | Disposal   | 16        |
| <b>6</b>  | <b>Dimensional drawings</b>                      | <b>17</b> |
| <b>7</b>  | <b>Accessories</b>                               | <b>18</b> |

# 1 Product description

## 1.1 Introduction

Thank you for buying this InPro® 7000-VP Series 2-electrode conductivity sensor from METTLER TOLEDO. InPro 7000-VP Series 2-electrode conductivity sensors are intended to be used to measure conductivity in low (high purity water) to medium conductivity water. The sensors are available in a wide selection of process connections to meet every application need.

### Markings

The printed markings on each InPro 7000-VP Series 2-electrode conductivity sensor contain the following information:

|                                 |  |
|---------------------------------|--|
| Mettler Toledo                  | Sensor manufacturer  |
| InPro 700X-VP                   | Type of sensor   |
| Cell M: XXXXXX                  | Cell constant  |
| Temp M: XXXXXX                  | Temperature constant   |
| Serial No: XXXXXXXXX            | Serial number  |
| Ex-classification               | Ⓔ II I/2 G Ex ia IIC T3/T4/T5/T6 Ga/Gb<br>(metallic sensors/metallic sensors with plastic surface) |
| Examination certificate number: | SEV 14 ATEX 0129 X   |
| Notified body:                  | Electrosuisse SEV, (No.: 1258)   |

## 1.2 Equipment and scope of delivery

### Scope of delivery

The InPro 7000-VP Series sensors are delivered ready for use. Each sensor is accompanied by this instruction manual, an individual «Certificate of Quality», CE Declaration of Conformity Certificate and a Material Certificate following EN 10204 (exceptions: InPro 7000-VP and InPro 7005-VP). Other Certificates may be included as specified for the individual product. Please check that the details given in the Certificate of Quality match the sensor label. For each sensor, the relevant cell constant has been determined individually in an ultrapure water system at 25 °C (77 °F) during the

manufacturing process, and the value documented in the accompanying «Certificate of Quality». It is not necessary to re-calibrate the sensor before initial operation.

### **Packaging**

The packaging consists of cardboard and plastic material.

Keep the packaging for later use during storage or transportation of the sensor. Should you wish to dispose of the packaging material, please observe your local regulations as well as those data and instructions given in Chapter 5.3 of this manual.

### **Unpacking and inspection**

Please check the sensors immediately during unpacking in order to determine possible damage or missing items. Any irregularities should immediately be reported to your carrier and to your supplier.

## 1.3 Technical data

|   | InPro 7000-VP                     | InPro 7005-VP                     | InPro 7001/120-VP<br>InPro 7001/225-VP            | InPro 7002/*-VP  |
|---|-----------------------------------|-----------------------------------|---|--|
| Measurement principle                                     | 2-electrode sensor                | 2-electrode sensor                | 2-electrode sensor                                | 2-electrode sensor   |
| Electrode material  | Titanium                          | Titanium                          | SS 316L   | SS 316L  |
| Body material   | PVDF                              | PTFE-coated<br>SS 316/1.4401      | SS 316L   | SS 316L  |
| RTD   | Built-in Pt 1000                  | Built-in Pt 1000                  | Built-in Pt 1000                                  | Built-in Pt 1000   |
| Sensor diameter   | See drawing                       | See drawing                       | 12 mm   | See drawing  |
| Insertion length  | 29 mm<br>(1.15")                  | 34 mm<br>(1.35")                  | 120/225 mm<br>(4.71/8.86")                        | 85/104 mm<br>(3.35/4.09")  |
| Max. sensor length  | 153.20 mm<br>(6.03")              | 75 mm<br>(2.95")                  | 194/299 mm<br>(7.64/11.77")                       | 156/175 mm<br>(6.14/6.88")   |
| Process connection  | - ¾" NPT<br>- 1" NPT Conduit      | - ¾" NPT                          | - Pg 13.5   | - Tri-Clamp 1.5"<br>- Tri-Clamp 2"<br>- Tuchenhagen-<br>Varivent<br>DN 40–DN 125 |
| Measuring range   | See separate table on page 8      |                                   |   |  |
| Cell constant nominal                                     | 0.1 cm <sup>-1</sup>              | 0.1 cm <sup>-1</sup>              | 0.1 cm <sup>-1</sup>                              | 0.1 cm <sup>-1</sup>   |
| Cell constant accuracy                                    | ± 1.0%                            | ± 1.0%                            | ± 1.0%  | ± 1.0%   |
| System accuracy <sup>a</sup>                              | ± 3.0% or better                  | ± 3.0% or better                  | ± 3.0% or better                                  | ± 3.0% or better   |
| <b>Working conditions</b>                                 |                                   |                                   |   |  |
| Max. pressure<br>at 25 °C (77 °F)                         | 34 bar<br>(500 psi)               | 17 bar<br>(250 psi)               | 17 bar<br>(250 psi)                               | 31 bar<br>(449.5 psi)  |
| Max. pressure<br>at 95 °C (203 °F) <sup>b</sup>           | 7 bar<br>(100 psi)                | 7 bar<br>(100 psi)                | 7 bar<br>(100 psi)                                | 10 bar<br>(145 psi)  |
| Measuring<br>temperature range                            | -10 ... 100 °C<br>(14 ... 212 °F) | -10 ... 100 °C<br>(14 ... 212 °F) | -10 ... 100 °C<br>(14 ... 212 °F)                 | -10 ... 120 °C<br>(14 ... 248 °F)  |
| Temperature range<br>(sterilization)                      | N/A                               | N/A                               | Sterilizable<br>-10 ... 131 °C<br>(14 ... 268 °F) | Sterilizable<br>-10 ... 155 °C<br>(14 ... 311 °F)                                |
| Temperature accuracy<br>at 25 °C (77 °F)                  | ± 0.25 °C<br>(± 0.5 °F)           | ± 0.25 °C<br>(± 0.5 °F)           | ± 0.25 °C<br>(± 0.5 °F)                           | ± 0.25 °C<br>(± 0.5 °F)  |
| <b>Design</b>   |                                   |                                   |   |  |
| Temperature<br>compensation                               | Pt 1000<br>IEC class A            | Pt 1000<br>IEC class A            | Pt 1000<br>IEC class A                            | Pt 1000<br>IEC class A   |
| Cable connection  | VarioPin (IP68)                   | VarioPin (IP68) <sup>c</sup>      | VarioPin (IP68)                                   | VarioPin (IP68)  |
| Wetted parts:   |                                   |                                   |   |  |
| - Metals  | Titanium (Grade 2)                | Titanium (Grade 2)                | SS 316L   | SS 316L  |
| - Plastics  | PVDF (FDA)                        | PTFE-coated<br>SS 316/1.4401      | -   | -  |
| - O-rings   | Viton® (FDA)                      | Viton® (FDA)                      | Viton® (FDA)                                      | Viton® (FDA)   |
| - Insulation  | PEEK (FDA)                        | PEEK (FDA)                        | PEEK (FDA)  | PEEK (FDA)   |
| - Surface roughness of<br>wetted metal parts <sup>d</sup> | N/A                               | N/A                               | Polished<br>Ra < 0.2 µm (< 8 µin)                 | Electropolished<br>Ra < 0.2 µm (< 8 µin)   |
| <b>Certificates</b>                                       |                                   |                                   |   |  |
| Cell constant   | •                                 | •                                 | •   | •  |
| Material certificate                                      |                                   |                                   |   |  |
| EN 10204 3.1  | -                                 | -                                 | •   | •  |
| Material specification 2.1                                | •                                 | •                                 | •   | •  |
| Surface roughness   | -                                 | -                                 | •   | •  |
| ATEX  | •                                 | •                                 | •   | •  |

### Captions

<sup>a</sup> The system accuracy depends on the measuring range and the particular transmitter.

<sup>b</sup> According to PED (Pressure Equipment Directive, 97/23/EC, Art. 3, Sec. 3).

<sup>c</sup> The VP is at the end of an approx. 0.5 m long fixed cable.

<sup>d</sup> Except at active electrode areas.

For each sensor, the cell constant and temperature constant have been measure/established individually and the values documented in the accompanying Quality Certificate. All calibrations can be traced back to NIST and/or ASTM Standards.

Both the measuring range and the system measurement accuracy strongly depend upon which type of transmitter is employed.

### Practical measuring range

| Transmitter            | InPro 7000-VP/<br>7005-VP | InPro 7001-VP | InPro 7002-VP |
|------------------------|---------------------------|---------------|---------------|
| <b>M300</b>            | 0.02–2,000                | 0.02–500      | 0.02–2,000    |
| <b>M300 ISM</b>        | –                         | –             | –             |
| <b>M400 4-w/2-w</b>    | 0.02–2,000                | 0.02–500      | 0.02–2,000    |
| <b>M700</b>            | 0.02–10,000               | 0.02–500      | 0.02–2,000    |
| <b>M800 1-channel</b>  | 0.02–2,000                | 0.02–500      | 0.02–2,000    |
| <b>System accuracy</b> | ±3 %                      | ±3 %          | ±3 %          |

## 2 Safety

### 2.1 Application compatibility

The wetted material parts of the sensor (several different materials come into contact with the sample medium) can under some circumstance be incompatible with the particular composition of the process medium and/or of the operating conditions. Responsibility to verify application compatibility lies wholly with the user.

The compatibility of different types of material are outlined on

**<http://www.coleparmer.com/techinfo/chemcomp.asp>**.

Mettler-Toledo GmbH, Process Analytics accepts no responsibility whatsoever for the correctness or accuracy of such details.

## 2.2 Proper utilization

METTLER TOLEDO InPro 7000-VP Series sensors are intended solely for the precise measurement of conductivity in aqueous and some limited non-aqueous solutions in industrial applications.

Any other use, or any operation over and above that intended by the manufacturer, are not recommended and can lead to harm or injury to material/equipment and persons. This is also relevant for applications which do not comply with the technical data labeled on the sensor. For any damage arising from such misuse, the user assumes full and sole responsibility.

## 2.3 Safety measures

The InPro 7000-VP Series sensors have been manufactured in line with state-of-the-art technology and in accordance with accepted technical safety regulations. Nevertheless, the sensors can still represent a source of risk and danger:

- if the sensors are operated by insufficiently trained personnel,
- if the sensors are employed incorrectly or not as intended by the manufacturer,
- if the sensors are not regularly maintained or serviced.

Local legislation and regulations must be observed at all times. Such stipulations do not form an integral part of this instruction manual.

It is necessary to use protective gear, including gloves, for persons coming in contact with the fluid. It is recommended to use gloves when handling the electrodes of the sensor to limit contamination of the sensor.

The user is responsible for the instruction and training of personnel. Additional copies of the instruction manual can be ordered from your supplier. This instruction manual is an essential element of the sensor equipment and must at all times be readily available to operators directly at the measurement site.

Before the sensor is removed from the process/process adapter, it must be ensured that the process pressure has been reduced to a safe level and the process temperature lowered to a safe range. Any escape of hot process fluid under pressure can cause damage to material/equipment or injury to persons.

No modification whatsoever may be carried out on the sensors. Any unauthorized modification or manipulation of the sensors results in immediate expiration of the full scope of warranty granted by the manufacturer.

## 2.4 Use in Ex-zones

### Note!

For intended installation in an Ex-classified area, please observe the following guidelines:

### Ex- classification:

II 1/2 G Ex ia IIB T3/T4/T5/T6 Ga/Gb (plastic sensors)

II 1/2 G Ex ia IIC T3/T4/T5/T6 Ga/Gb (metallic sensors/metallic sensors with plastic surface)

### Marking and number of the test certificate:

SEV 14 ATEX 0129 X

### Introduction

1. The conductivity sensors made of plastic are according to Directive 94/9/EC (ATEX 95) Annex I equipment group II, category 1/2G after RL 99/92/EC (ATEX 137) in the zones 0/1 or 0/2 or 1/2 and the gas groups IIA and IIB, explosion hazards due to flammable substances in the field of temperature classes T1 to T6, may be used.
2. The conductivity sensors made of metal or metal coated with plastic to Directive 94/9/EC (ATEX 95) Annex I to equipment group II, category 1/2G after RL 99/92/EC (ATEX 137) in zones 0/1 or 0/2 or 1/2 and gas groups IIA, IIB and IIC, explosion hazards due to flammable substances in the field of temperature classes T1 to T6, may be used.
3. For use/installation, the requirements according to EN 60079-14 must be observed.

## Electrical ratings

Conductivity measuring circuit and temperature sensing circuit in type of protection intrinsic safety Ex ia IIC only for connection to a certified intrinsically safe circuit.

Maximum values:

|       |        |        |
|-------|--------|--------|
| $U_i$ | $\leq$ | 16 V   |
| $I_i$ | $\leq$ | 190 mA |
| $P_i$ | $\leq$ | 200 mW |

The effective internal inductance and capacitance are negligible.

The values above apply, each as the sum of all the individual circuits of the associated intrinsically safe supply and evaluation unit.

## Special conditions for safe use

1. The relationship between the maximum permissible ambient or media temperature and temperature class is shown in the following table:

For  $U_i \leq 16$  V,  $I_i \leq 190$  mA,  $P_i \leq 200$  mW.

Conductivity measuring circuit and temperature measuring circuits:

| Temperature class | Maximum ambient or media temperature |
|-------------------|--------------------------------------|
| <b>T6</b>         | 51 °C                                |
| <b>T5</b>         | 63 °C                                |
| <b>T4</b>         | 91 °C                                |
| <b>T3</b>         | 143 °C                               |


2. For installation, the system drawings and wiring diagrams of the measurement systems should be observed.
3. The capacitance and inductance of the connecting cable has to be considered.

4. The metal body of the conductivity sensors or the safety weld-in socket or the fitting InFit 76\*\_\*\_\*\_\* or InTrac 7\*\*\_\*\_\*\_\* or other appropriate fitting must be electrically connected to the equipotential bonding system.
5. The metal body of the conductivity sensors or the safety weld-in socket or the fitting InFit 76\*\_\*\_\*\_\* or InTrac 7\*\*\_\*\_\*\_\* or other appropriate fitting is optionally included in the routine pressure test of the system.

## 3 Initial start-up

### 3.1 Installation

The sensors are to be mounted in such a way that the sample medium flows directly into the sensor through the orifice at the tip and re-emerges through the vent holes. Any other installation position of the sensors can create the risk of formation of airlocks, or of contamination through deposits of solid matter.

Flow of sample medium must be toward the face of the sensor. Air bubbles and solid matter deposits are to be avoided. A minimum clearance of 0.6 cm (1/4") between sensor and pipe wall must be maintained at the tip of the sensor.

Vertical mounting (top entry) in a pipe is recommended only if the pipe is full and no air bubbles are able to develop in the flow.

In the event of side mounting (side entry) of the sensor, vertical upward flow of sample process medium must be ensured.

## 3.2 Electrical connections

All of the InPro 7000-VP Series sensors can be connected to the associated transmitter using the appropriate multiwire cable from METTLER TOLEDO. Depending on which type of transmitter is employed, not all cable strands may have to be used.

### Colors and function of the cable strands (supplied separately)

Connection to the relative transmitter must be followed according to the wiring diagram below.

| <b>2-electrode sensor</b> | <b>Color</b> |
|---------------------------|--------------|
| Inner electrode           | white        |
| Outer electrode           | white/blue   |
| RTD                       | green        |
| RTD sense                 | red          |
| RTD/GND                   | Bare shield  |

Connection to the relative transmitter must be followed according to the wiring diagram (see respective transmitter instruction manual for detail).

METTLER TOLEDO supplies the cables. Please refer to chapter 7 "Accessories" for ordering information

## 4. Maintenance and troubleshooting

### 4.1 Conditions of warranty

METTLER TOLEDO guarantees the quality of materials and workmanship within a narrow range of manufacturing tolerances, so that the product purchased is free from any substantial deviations from material and manufacturing quality standards. The warranty is valid for a period of one year from date of delivery: If within this warranty period, any repair or replacement should become necessary, and such cause is not due to misuse or incorrect application, please return the sensor, freight pre-paid, to your appropriate METTLER TOLEDO supplier. Repair work will be carried out free of charge. Final decision on whether the defect is due to a manufacturing error or to incorrect operation of the sensor by the customer is made at the option of the Customer Service department of METTLER TOLEDO. After expiration of the warranty period, sensors will be repaired or replaced on an exchange basis against payment of the costs involved.

### 4.2 Maintenance

Dirty or contaminated sensors can deliver incorrect measurement values. If fouling is presumed, the sensor is to be removed from operation and the electrodes as well as the insulation between the electrodes cleaned with a soft cloth. Suitable cleaning solutions are mild detergents or diluted acids (< 0.5 % by wt.) such as nitric acid. Never use cleaning agents that are not compatible with the material to be cleaned.

When handling acids, precautionary measures are to be taken.

The sensor must be thoroughly flushed with distilled or deionized water prior to re-installation.

Following cleaning and re-installation into the process, it can take from several minutes up to several hours until the sensor delivers the originally measured value.

## 4.3 Troubleshooting

| <b>Error</b>                              | <b>Possible cause</b>   | <b>Corrective action</b>   |
|---|---|--|
| No signal on display                      | Electrical connections either missing or incorrectly in place | Check all connections and associated cabling   |
|   | Sensor is not in contact with the sample medium               | Check installation for air pockets   |
| No temperature signal                     | Temperature probe is not connected                            | Check all connections and associated cabling   |
|   | Transmitter does not support Pt1000 RTD                       | Replacement of transmitter or sensor necessary   |
| Incorrect or unstable measurement reading | Sensor is contaminated/fouled                                 | Clean the sensor (electrodes)  |
|   | Sensor is installed too close to pipe wall                    | Ensure minimum distance of 0.6 cm (1/4") between sensor tip and pipe wall is maintained. |

## **5 Removal from operation, storage, disposal**

### **5.1 Removal from operation**

The sensor is only conditionally subject to aging. When used as intended and appropriately maintained and serviced, the lifetime of the sensor can extend to several years.

Before removing the sensor from the process/process adapter, it must be ensured that the process pressure has been reduced to a safe level and the process temperature has lowered to a safe range. Any escape of hot process medium under pressure can cause damage to material/equipment or injury to persons.

After removal from the process the sensor should first be flushed with distilled water.

If the sensor incurs a defect, it can not be repaired and must be disposed under observance of prevailing local regulations.

### **5.2 Storage**

If the sensor is not use, it can be stored dry. However, it has to be reconditioned in the process fluid accordingly prior to renewed operation. This procedure may take several hours to complete.


### **5.3 Disposal**

The user is responsible for the proper disposal of the conductivity sensor. The sensor contains electronic components that must be disposed of correctly in order to avoid hazards to people or to the environment.


## 6 Dimensional drawings


**InPro 7000-VP**


**InPro 7005-VP**


**InPro 7001/120-VP**  
**InPro 7001/225-VP**

**InPro 7002/1.5"-TC-VP**  
**InPro 7002/2"-TC-VP**

**InPro 7002-VAR-VP**

## 7 Accessories

### O-rings

| Sensor                                 | Standard O-ring             | Complementary O-rings   | Order No.                              |
|--|-----------------------------|---|--|
| InPro 7001/120-VP<br>InPro 7001/225-VP | Viton FDA<br>(10.77 x 2.62) |   | 59 905 680                             |
|  |                             | Silicon FDA (10.77 x 2.62)<br>Kalrez USP VI (10.77 x 2.62)<br>EPDM FDA (10.77 x 2.62) | 59 905 678<br>59 905 789<br>59 905 768 |

### Cable

The max. cable length recommended is 60 m.

| Cable length    | Order No.  |
|-----------------|------------|
| 1.5 m (5 ft)    | 58 080 201 |
| 3.0 m (10 ft)   | 58 080 202 |
| 4.6 m (15 ft)   | 58 080 203 |
| 7.6 m (25 ft)   | 58 080 204 |
| 15.2 m (50 ft)  | 58 080 205 |
| 22.9 m (75 ft)  | 58 080 206 |
| 30.5 m (100 ft) | 58 080 207 |

### Housings for InPro 7000-VP Series

| Designation | Applied as          | Suitable sensor             |
|-------------|---------------------|-----------------------------|
| InFit 761 e | Stationary housing  | InPro 7001/120 and 225-VP   |
| InFit 762 e | Stationary housing  | InPro 7001/120-VP           |
| InFit 777 e | Retractable housing | InPro 7001/225-VP           |
| InFit 787 e | Retractable housing | InPro 7001/120-VP           |
| InFit 798 e | Retractable housing | InPro 7001/120-VP           |
| InDip 550   | Immersion housing   | InPro 7000-VP/InPro 7005-VP |

Please contact your METTLER TOLEDO representative for more options.

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**France Mettler-Toledo Analyse Industrielle S.A.S.**
**Germany Mettler-Toledo GmbH**
**Switzerland Mettler-Toledo (Schweiz) GmbH**