PREGNANCY DETECTOR (PD)

MANUAL

EN

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DE

MODE D'EMPLOI

FR

MANUAL DE USO

ES

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CE

TABLE OF CONTENTS

INTRODUCTION	3
DESIGN	5
CHECKING THE DEVICE AND FIRST MEASUREMENTS	7
MEASUREMENTS IN ANIMALS	11
FINAL NOTES FOR USERS	17
BATTERY REPLACEMENT	19
TECHNICAL DATA	21

EN

INTRODUCTION

EN

1

CHAPTER

The tester is meant for breeders, veterinarians, inseminators and advisory services experts. It is very helpful in situations when it's necessary to detect pregnancy fast in its early phase. ULTRASOUND PREGNANCY DETECTOR is used to confirm pregnancy in sows, goats and ewes by detecting amniotic fluid in the uterus.

The probe applied to the animal's skin sends an ultrasound beam which is reflected from uterine horn filled with fluid and returns to the probe, causing a change of signal generated by the device. When pregnancy is detected, the sound starts to pulsate more intensely and a green diode is flashing.


ULTRASOUND PREGNANCY TESTER was designed in such a way to operate it easily and perform fast and effective examinations. The device is completely safe for animals' health.

DESIGN

EN

2

CHAPTER

The device is small, light and handy, which makes it easy to use.

CHECKING THE DEVICE AND FIRST MEASUREMENTS

EN

3 **CHAPTER**

A red diode signals that the battery has discharged.

A green diode is closely linked with the tone signalling. When the power is on, the diode will shine.

When the front of the ultrasound head comes into a good contact with the animal's body, the diode will begin to flicker with the same frequency as the tone signal – about one impulse per second.

Upon detecting pregnancy, diode flickering and frequency of the tone emitted by the sound generator will increase significantly - up to 4-5 impulses per second.

The ultrasound probe is small and handy, so that it can be manipulated easily.

Preparing the detector

At the beginning of each examination, a cable connecting the probe with the device should be carefully connected. Next, the power should be turned on - this will cause the diode to shine a continuous light, which indicates that the device is operable and ready to use.

If during the operation or directly afterwards, the red diode starts to shine, in addition to the green diode, this means that the battery should be replaced.

The device has an automatic switch off system. After 3-4 minutes from switching on the device, it will switch off automatically. It allows for saving the battery in case of accidental start-up or forgetting to switch off the device between or after examinations.

A test to check the operation of the device

After connecting the probe and switching on the device, a test to check its operation may be performed:

- Moisten the head of the probe and then apply it to the inner side of your forearm, several centimetres from the wrist.
- Manipulate the probe as long as an intermittent signal is heard, which indicated good contact with the skin (pulsating light and sound - about 1 impulse / sec.). Insufficient amount of water may cause lack of such signal.
- Put the probe tip to a glass of water (at least half a glass), pointing the ultrasound beam towards the bottom - the device will start up the same signal as when pregnancy is detected (intensely pulsating light and sound, 4-5 impulses per second).

Any other reaction of the device may indicate damage. In such a case you must contact the Dramiński service centre.

MEASUREMENTS IN ANIMALS

EN

4 **CHAPTER**

Before each examination, date of mating should be checked up in the documentation and a likely pregnancy date specified, which will help to perform the test at the right time.

- The tip (head) of the probe should be immersed in moisturising liquid (vaseline oil, liquid paraffin, cooking oil or technical oil).
- Moisture the animal's skin in the appropriate place. It is necessary for the probe to have a good contact with the animal's skin.
- Try to perform the examination at the right time (see the table).
- If the examination is performed at a very early stage of pregnancy, it definitely must be repeated at the right time (see the table).

Animal	Pregnancy detector	
	Date of detection (possible)	Date of detection (optimum)
Sow 	18-75	25-35
Ewe 	40-120	60-90
Goat 	40-120	60-90

How to press the probe


- The probe should be held quite firmly, however it shouldn't be pushed too hard in order to avoid significant skin deformation.
- Point the ultrasound beam towards the head at the right angle (see the picture). Pointing the beam close to a vertical position or backwards may cause the beam to reflect from the full bladder and generate a false pregnancy result.
- If the moistened probe is placed correctly, the device will signal good contact with the animal's body by changing the light of the green diode from continuous into regularly flickering one (frequency of 1 second) and a tone signal pulsating with the same frequency.
- This is when you should delicately move the probe in different directions to penetrate the abdominal

cavity, where ultrasounds will meet the uterine horn filled with amniotic fluid.


- If amniotic fluid is detected, the light and sound signal will start pulsating intensely (with the frequency of 4-5 times per second).

If the result of the test in a particular place is negative (the signal indicating good contact fails to change into a more intense one), the examination should be repeated in another place, first closer to the rear leg and then further from it (remembering that both the probe and the skin need to be well moistened).

If the pregnancy signal is not obtained despite careful penetration of the abdominal cavity, the result of the examination should be considered negative, provided that it has been such for three consecutive days.

Ewe

Pig

The examination should be performed on a standing animal, at both sides of its body, starting from the right hand side.

The examination takes up more or less a dozen seconds, or several seconds after some practice.

As far as possible, the animal should be calm and relaxed during the examination.

The best time to perform the examination is in the morning, before the animals have eaten and after they have urinated. Full bladder and stomach (watering the animal before the examination) may lead to obtaining a false result (ultrasound beam being reflected not from amniotic fluid but from full bladder or stomach).

Repeating the test for three consecutive days ensures that the diagnosis is reliable.

FINAL NOTES FOR USERS

EN

CHAPTER

5

- You should avoid moistening and wetting the cable connector and the probe snap-ons.
- The device may be cleaned with a delicately moistened cloth, but it must be completely dry before switching it on.
- The device is designed to make the examination procedure easy.
- The more examinations are performed, the greater the skills of the user become, and pregnancy detection becomes fast and effective.

BATTERY REPLACEMENT

EN

6

CHAPTER

The red diode signals that the battery has discharged. It indicates that examinations should be discontinued immediately (false signals are possible) and the battery should be replaced before conducting further research.

At the back of the device, there is a flap that needs to be loosened. Remove the old battery and replace it with a new one. Next, fix the flap, paying attention to the gasket.

TECHNICAL DATA

EN

9

CHAPTER

Dimensions	Height - 4 cm Width - 8 cm Length - 15,5 cm Length of the cable – 110 cm
Weight of the device	330 g
Power	9V type 6LF22 alkaline battery
Operating and storing temperature	From 0° to + 60° C
Type of casing	ABS material
Fittings	Ultrasound pregnancy detector A probe with a connecting cable 9V alkaline battery A bottle of fluid paraffin Transport packaging Manual with guarantee

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